

Dissolved gases account for many of the common problems encountered in routine analyses, such as bubble formation, pump cavitation, detector noise, baseline drift, and gradient precision. The Finnigan SpectraSYSTEM solvent degassers ensure the optimum performance of your HPLC system by thoroughly removing dissolved gases from the mobile phase.

Finnigan™ SpectraSYSTEM™ Solvent Degassers

Ensure optimal performance
of your HPLC pumping system



- **More stable baselines, improved accuracy and precision**
- **Rapid degassing of up to four solvent channels**
- **Choice of vacuum degassing, helium sparging, or both**
- **Manual, timed, or programmed operation**
- **Biocompatible materials protect sensitive samples**
- **Leak sensor, automatic shutdown, and double-wall containment enhance lab safety**
- **Optimum sensitivity and selectivity with Thermo Electron columns**

Optimize Analytical Data

Solvent degassing gives you more stable baselines. Expensive HPLC-grade solvents minimize absorbance by the mobile phase. But, dissolved gases (e.g. oxygen) increase both short-term noise and long-term drift. Outgassing of dissolved gases in the flowcell causes detector spikes. Vacuum degassing is thus an effective, cost-efficient way of improving the quality of your HPLC data.

Improve both gradient shape and reproducibility. The compressibility of dissolved gases creates small volumetric inaccuracies which affect gradient shape and precision. Removal of these gases improves flow accuracy and precision for maximum pump performance.

Rapid Degassing of up to Four Solvent Channels

Thin-wall membrane tubing and high-capacity vacuum chambers permit fast, efficient removal of dissolved gases. For high flowrate applications, like semi-prep scale, a "T" connector is used to combine the output of individual solvent channels.

Match a Degassing Technique to Your Application

Thermo Electron offers customers helium sparging or vacuum degassing, or both. The high-efficiency vacuum degassing system suffices for the vast majority of applications. Choose the technique best suited to your application needs and budget, or choose both for maximum flexibility and performance.

The vacuum degasser is mounted in an attractive and ergonomic, vertical package. For helium sparging, you can mount up to a 4-channel SHM helium degassing kit on a vertical rail attached to the degasser box or onto your Finnigan SpectraSYSTEM pump.

Continuous Operation

The vacuum degasser operates continuously to ensure a stable level of vacuum through a sequence of analyses. This provides reproducible gradient or isocratic run-to-run performance.

Biocompatible Protection for Sensitive Samples

All wetted materials in the degasser are chemically and biologically inert. This ensures maximum corrosion resistance and compatibility with sensitive biomolecules.

Built-in Safety Features

A leak sensor and automatic shutdown enhance lab safety. The leak detection system gives an immediate warning of a leak. A red warning light activates, and the safety interlock turns off all power.

Double-wall containment enhances lab safety. In some environments, double-wall containment is required when the chance of a hazardous material leak exists. This unit satisfies all known safety requirements.

Vacuum Degasser

Method

Vacuum Membrane

Solvent Channels

Four

Solvent Contact Materials

Teflon® AF, TFE and titanium

Volume Per Channel

500 µL

Free Flow Per Channel

5 mL/min H₂O with 2 ft. head pressure

Dimensions

36 × 15 × 55 cm (H × W × D)

Weight

13 kg

Power Requirements

100, 120, 220, 240 VAC; 50/60 Hz

Product Certification

TUV CUS

TUV CE Mark

FCC (EMI)

Part Number Description

SCM1000-011	SCM1000 115V
SCM1000-022	SCM1000 230V
SCM1000E-011	SCM1000 Degasser 4L Bottle
A4122-010	Helium Degas Manifold (2)
A4144-010	Helium Degas Manifold (4)
A4125-010	Helium Degas Mounting Bracket
A5765-010	SpectraSYSTEM Accessory Rack

International Offices

Australia

Tel. (61) 2 9898 9000

Austria

Tel. (43) 1 333 50340

Belgium

Tel. (32) 2 482 3030

Canada

Tel. (905) 712 2258

Canada only

Tel. (800) 721 4260

Central and South America

Tel. (01) 512 251 1530

France

Tel. (33) 1 60 92 48 00

Germany

Tel. (49) 6103 4080

Italy

Tel. (39) 02 950 591

Japan

Tel. (81) 45 453 9100

The Netherlands

Tel. (31) 76 5878 722

People's Republic of China

Tel. (86) 10 6621 0839

Spain

Tel. (34) 91 657 4930

Sweden

Tel. (46) 8 556 468 00

Switzerland

Tel. (41) 61487 8400

United Kingdom

Tel. (44) 1 442 233555

United States

Tel. (01) 800 532 4752