Reference Emulsion 100240

I. Intended Use

Pacific Hemostasis® Reference Emulsion is intended to be used as a priming solution on Instrumentation Laboratory ACL™ coagulation analyzers and to provide an optical reference for coagulation testing.

II. Summary and Principles

The Instrumentation Laboratory ACL is a coagulation analyzer that detects the formation of a fibrin clot by monitoring the change in light scatter as a function of time. The Reference Emulsion functions as an optical blank that must be present while the instrument is performing optical reading for each test. The ACL instrument will automatically place a fixed amount of Reference Emulsion into one cuvette well and use it to normalize the light scatter from the other cuvette wells which contain patient samples. The ACL instrument will not report a result without the optical reading from the Reference Emulsion solution. ¹

III. Reagent

For *in vitro* diagnostic use

Composition:

Contains silicone oil emulsion (0.1%), antifoam (0.02%) and sodium azide (0.1%) in water. The turbidity of each lot has been carefully adjusted to provide optimum performance.

Warning. Reference Emulsion contains sodium azide. Sodium azide under acidic conditions yields hydrozoic acid, an extremely toxic compound. Azide compounds should be diluted with running water before being discarded. Upon disposal, azide compounds should be flushed with large amounts of water. These precautions are recommended to avoid deposits in metal pipes in which explosive condition may develop.

Storage Instructions and Stability

Store at 15-30 °C. **DO NOT FREEZE.**

Unopened containers of Reference Emulsion are stable until the expiration date printed on the product label.

IV. Specimen Collection

Refer to the instructions provided by the reagent and instrument manufacturers.

V. Test Procedure

Materials Required But Not Provided

The materials required but not provided include: ACL Coagulation analyzer, Coagulation reagents, and Coagulation controls.

Invert bottle 3-4 times before placing on the instrument. Follow the instrument manufacturer's directions for priming the instrument.

PT-FIB Menu Option. If the fibrinogen (FIB) concentration is being determined as part of the PT (prothrombin time) and the lot number of Reference Emulsion is different from the previous lot, the instrument must be calibrated. Follow the instrument manufacturer's directions and allow the instrument to calibrate the PT/FIB. Calibration is not required for PT or APTT.

VI. Limitations

The Reference Emulsion is intended to be used on only centrifugal-based ACL instrumentation.




All samples, reagents and controls are subject to the limitations of the test system. Variables such as temperature, reagent stability, instrument performance, and individual technique can influence final results. Always follow instrument and reagent manufacturer guidelines. The reagents and controls described here are the Pacific Hemostasis brand. The performance of other brands of reagents and controls may differ slightly and their performance must be validated by the user.

VII. Performance Characteristics

The following illustrates the expected precision when using Reference Emulsion with representative Pacific Hemostasis Reagents and controls.

Test	Control	Within Run (%)	Day-to-Day (%)
PT (Thromboplastin-DL)	Normal	1.3	0.8
	Abnormal	2.1	1.7
APTT (APTT-XL)	Normal	1.0	1.4
	Abnormal	0.6	4.6
Fibrinogen (derived from the PT)	Normal	2.7	1.4
	Abnormal	3.8	2.7

The following is a comparison of PT, APTT and Fibrinogen results when multiple samples are tested on an ACL 100 instrument primed with Reference Emulsion or a different brand.

VIII. References

- Operator's Manual for the ACL 3000

Ordering Information		
Cat. No.	Description	Contents
100240	Reference Emulsion	4 x 500 mL
100355	Thromboplastin-DL, 200 Det.	10 x 4 mL
100359	Thromboplastin-DL, 500 Det.	10 x 10 mL
100402	APTT-XL, 400 Det.	10 x 4 mL
100403	APTT-XL, 1000 Det.	10 x 10 mL
100304	Calcium Chloride- 0.020 M	10 x 10 mL
100595	Coagulation Control, Level 1	10 x 1 mL
100596	Coagulation Control, Level 2	10 x 1 mL
100597	Coagulation Control, Level 3	10 x 1 mL
100604	Low Fibrinogen Control	10 x 1 mL
100607	High Fibrinogen Control	10 x 1 mL

Fisher Diagnostics® Limited Warranty

Fisher Diagnostics® (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

Pacific Hemostasis® is a registered trademark of Fisher Scientific Company L.L.C.

Fisher Diagnostics® is a registered trademark of Fisher Scientific Company L.L.C.

ACL™ is a registered trademark of Instrumentation Laboratory.

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MDCI Ltd.