## Immunodepleted Factor VIII Deficient Plasma (Human)

### I. Intended Use

Pacific Hemostasis® Immunodepleted Factor VIII Deficient Plasma is intended for use as a substrate in the quantitative determination of Factor VIII activity in citrated plasma.

### II. Summary and Principles

Factor VIII activity in plasma is assayed by the amount of Activated Partial Thromboplastin Time (APTT) correction produced by the test plasma, when mixed with factor deficient plasma. The correction of the unknown is compared to that produced by a reference plasma of known normal factor activity, such as Pacific Hemostasis Universal Coagulation Reference Plasma (UCRP™).

### III. Reagent

For *in vitro* diagnostic use.

#### Composition:

A lyophilized plasma containing stabilizer depleted of Factor VIII by immobilized highly specific antibodies. Factor VIII activity is less than 1% of normal levels; all other coagulation factors are within the normal range. Store unopened vials at 2-8°C. Each reagent is stable at 2-8°C until the expiration date printed on the vial label. Reconstitute with 1.0 mL of distilled water. Agitate gently until solution is complete. **Do not freeze reconstituted plasma.** The reconstituted plasma is stable for 8 hours when stored stoppered at 2-8°C. Erratic values, product color variations, or lack of vacuum in the vials could indicate product deterioration.

**CAUTION! POTENTIAL BIOHAZARD:** Contains human source material. While each human serum or plasma donor unit used in the manufacture of this product was tested by FDA-approved methods and found nonreactive for hepatitis B surface antigen (HBsAg), antibody to hepatitis C (HCV), and antibody to HIV-1/2, all products manufactured using human source material should be handled as potentially infectious. Because no test method can offer complete assurance that hepatitis B or C viruses, HIV, or other infectious agents are absent, these products should be handled according to established good laboratory practices.

### IV. Specimen Collection

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 2500 x g. Test within 2 hours if samples are held at 22-24°C. For more details on specimen collection and storage, see NCCLS Document H21-A3.1

### V. Test Procedure

This product is suitable for use with manual, mechanical, photo-optical or nephelometric methods for clot detection. Consult instrument manufacturer for more specific guidelines. The following procedural outline assumes use of a manual method.

#### Materials Provided:

Immunodepleted Factor VIII Deficient Plasma (Human, Dried), 10x1 mL

#### Materials Required but Not Provided:

- Pacific Hemostasis APTT reagent.
- Pacific Hemostasis Calcium Chloride reagent.
- Pacific Hemostasis Diluting Fluid (Barbital Buffered Saline).
- Pacific Hemostasis Universal Coagulation Reference Plasma.
- Plastic test tubes. Precision pipettes.
- Distilled water. Log-log graph paper.

#### Notes on Procedure:

- Perform all testing in duplicate. Duplicates should agree within 5%, repeat if necessary.
- Prepare dilutions immediately prior to testing.
- A new standard curve must be prepared each time testing is performed.

#### Preparation of Factor Activity Reference Curve:

- Reconstitute one vial of Pacific Hemostasis UCRP. See package insert accompanying this reagent for reconstitution instructions, storage, and stability.
- Reconstitute and pool sufficient vials of the Immunodepleted Factor Deficient Plasma to allow 0.2 mL for each dilution to be tested.
- To construct the standard curve, label 5 plastic test tubes and pipet Pacific Hemostasis Diluting Fluid (Barbital Buffered Saline) and UCRP into each tube according to the chart below. Cap and mix gently.

Tube#	Diluting Fluid	UCRP	Dilution	Dilution Factor	% Factor Activity in Dilution*
1	0.9 mL	0.1 mL	1:10	1	110
2	1.9 mL	0.1 mL	1:20	2	55
3	3.9 mL	0.1 mL	1:40	4	27.5
4	3.95 mL	0.05 mL	1:80	8	13.8
5	1.0 mL	1.0 mL tube #4	1:160	16	6.9

\*Chart is example only. Based on assigned value of 110% for UCRP factor activity.

#### Test Standards, Controls and Test Samples:

- Prepare dilutions of patient and control samples. At least three different dilutions of the patient samples should be tested, beginning with the same dilution as the first point of the curve (i.e. 1:10, 1:20, 1:40).
- Prewarm Pacific Hemostasis Calcium Chloride reagent to 37°C. Refer to package insert with APTT reagent for information on the appropriate Calcium Chloride to use.
  - Pipet 0.1 mL APTT reagent into a cuvette. Incubate at least 3 minutes, but no longer than 20 minutes.
  - Add 0.1 mL Factor Deficient Plasma and 0.1 mL of standard, control or patient plasma dilution to the cuvette. Mix well. Refer to APTT reagent package insert for length of activation.
  - Add 0.1 mL prewarmed Calcium Chloride and time clot formation.


### VI. Results

#### Reference Curve

- For each dilution of the UCRP, calculate the mean of the duplicate clotting times to the nearest 0.1 second.
- Calculate the factor activity present in each dilution of the UCRP. The first dilution (1:10 in the example) corresponds to 100% factor activity; therefore the factor activity of this dilution is equal to the assay value of the UCRP. The dilution factor represents the relationship between the first dilution and subsequent dilutions. Divide the assay value of the standard by the dilution factor to determine the percent activity for each dilution. Refer to the table in section V.
- Plot the mean clotting time of each dilution against its percent activity using log-log graph paper. Construct a straight line of best fit.

#### Test Specimens and Controls

- For each dilution of the test specimen, calculate the mean of the duplicate clotting times to the nearest 0.1 second.
- Locate the mean clotting time on the vertical axis of the reference curve. Find the corresponding point on the reference line and read the percent activity of the dilution on the horizontal axis of the graph. Clotting times outside the curve boundaries should not be used.
- Multiply the percent activity from the graph by the appropriate dilution factor to determine the actual percent activity for each dilution. If the values recovered for the individual dilutions do not agree with each other, new dilutions should be made and the testing repeated. If they still do not agree, you should consider the possibility of an inhibitor.

Sample Calculation. Using the above reference curve, a 1:10 dilution of patient plasma that yields 60.0 seconds would have a calculated factor activity of 10%. A 1:20 dilution that yields 60.0 seconds would have a factor activity of 20%. (Dilution factor x activity = final result).

### VII. Quality Control

Assayed reference plasmas, such as Pacific Hemostasis Abnormal Coagulation Reference Plasma (ACRP™) and an additional lot number of UCRP should be tested to validate the reference curve. ACRP is assigned values for coagulation factors in the abnormal range. Each lot of ACRP will vary, but the values are expected to be approximately 30-50% of normal. Use of this plasma provides an excellent quality control check on the accuracy of factor determinations at the low end of the reference curve. Testing an additional lot of UCRP (other than the one used for the standard curve) should be used as a control in the normal range. Actual values recovered depend on individual technique, instrument, standard, and reagent used.

### VIII. Limitations

Variables such as temperature, reagent stability, plasma sample conditions, instrument performance, as well as individual technique can influence factor activity results. Since there are many variables contributing to assay performance, each laboratory must establish their own reference range, as well as acceptable limits of quality control. Discrepant results between dilutions may be indicative of the presence of an inhibitor.

### IX. Expected Values

Normal plasma will yield factor activity values between 50-150% of normal.<sup>4,5</sup>

### X. Performance Characteristics

Pacific Hemostasis Immunodepleted Factor VIII Deficient Plasma will perform according to the limitations of the procedure described herein when tested with Pacific Hemostasis reagents. Precision and accuracy of the deficient plasma was assessed by preparing standard curves on 8-10 days, and determining the recovery of Factor VIII activity contained in five frozen control plasmas. The FVIII values recovered were within 2% of the expected control values. The results are shown below: <sup>1-3</sup>

	Day-to-Day Recovery of FVIII Activity (%) Contained in Control Plasmas*				
	Control #1	Control #2	Control #3	Control #4	Control #5
mean	130.9	95.5	65.6	33.6	18.0
1 SD	6.46	4.63	1.82	2.13	1.32
% CV	4.9	4.8	2.8	6.3	7.3
days of testing	8	8	10	8	8

\*Testing with Pacific Hemostasis APTT-LS reagent on the MLA®-1000C™. Standard curve R<sup>2</sup> = 0.99 on all days of testing.

### XI. References

- NCCLS. *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays*. 3rd edition, Approved guideline. NCCLS Document H21-A3. Wayne, PA, 1998.
- Powers, L.W. *Technical Hematology*. In: Diagnostic Hematology. 1989. Bircher S. (Ed.). The C.V. Mosby Co., St. Louis, MO, p.484.
- Stability data found in DHF.
- 4-5. Data found in 510K file.

Ordering Information		
Cat. No.	Description	Contents
100806	Immunodepleted FVIII Deficient Plasma	10 x 1 mL
100644	Diluting Fluid (Barbital Buffered Saline)	2 x 100 mL
100301	UCRP	10 x 1 mL
100603	ACRP	10 x 1 mL

### FISHER DIAGNOSTICS® LIMITED WARRANTY

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

Pacific Hemostasis® is a registered trademark of Fisher Scientific Company L.L.C. Fisher Diagnostics® is a registered trademark of Fisher Scientific Company L.L.C. MLA®-1000C™ is a registered trademark of Instrumentation Laboratory. UCRP™ is a registered trademark of Fisher Scientific Company L.L.C. ACRP™ is a registered trademark of Fisher Scientific Company L.L.C.

MDCI Ltd.  
Arundel House  
1 Liverpool Gardens  
Worthing, West Sussex BN11 1SL  
UK

Fisher Diagnostics®  
Fisher Scientific Company L.L.C.  
8365 Valley Pike  
Middletown, VA 22645-0307, USA  
1-(800) 528-0494 USA only

## Immunodepleted Factor VIII Deficient Plasma (Human)

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