## Low Fibrinogen Control Plasma or High Fibrinogen Control Plasma

### I. Intended Use

Pacific Hemostasis® Low and High Fibrinogen Control Plasmas are intended for use in the quality control of fibrinogen assays.

### II. Summary and Principles

The use of controls is a valuable tool for assuring quality in coagulation testing. Fibrinogen Controls contain a known amount of fibrinogen in the low and high abnormal range. These controls provide reliable markers for monitoring the accuracy of abnormal fibrinogen results and can be used to verify the extent of the reportable range.

### III. Reagent

For *in vitro* diagnostic use.

**Composition:** Fibrinogen Controls: Processed from human plasma collected with sodium citrate anticoagulant (4% w/v). The plasmas are adjusted to yield fibrinogen values in the low as well as high abnormal range. <1.0% Stabilizers and buffers are added prior to lyophilization.

Store unopened vials at 2-8°C. Reconstitute with 1.0 mL of distilled water. Swirl gently and let stand undisturbed for 15 minutes at room temperature. Do not invert vial or mix vigorously. After proper reconstitution, High Fibrinogen Control is stable for 8 hours when stored in capped vial at 2-8°C and Low Fibrinogen Control is stable for 16 hours when stored in capped vial at 2-8°C.<sup>3</sup> A precipitate may form when refrigerated. Gently warm the plasma to 37°C to minimize any precipitate.

Erratic values, product color variations, or lack of vacuum in the vials could indicate product deterioration. However, poor control performance could also be due to other factors within the test system.

**Caution:** Each unit of source material used in the preparation of these products has been tested by an FDA licensed method and found non-reactive for HBsAg and negative for antibodies to HIV and HCV. However, no known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS, or other infectious diseases. These products, like all materials of human origin, should be handled as potentially infectious biological material.

### IV. Test Procedure

**Materials Provided:** Low Fibrinogen, 10x1 mL or High Fibrinogen Control Plasma, 10x1 mL

**Materials Required, But Not Provided:**

Distilled water  
Pipette capable of accurately delivering 1.0 mL

Test Low and High Fibrinogen Controls following protocols established for patient samples. A 1:20 dilution is recommended for High Fibrinogen Control. Actual results depend on many factors, including lot number, reagent, and instrument used. Ranges must be determined in each laboratory with changes of lot number, reagent, or instrument. Normal and abnormal controls should be run at the initiation of testing each day and at least once each shift, or with

each group of assays. Controls should also be tested with each reagent change or major instrument adjustment. Follow instructions provided by instrument and reagent manufacturers.

### V. Results

Results depend on reagents, instruments, and individual laboratory protocols. However, Low Fibrinogen and High Fibrinogen Controls are formulated to have fibrinogen values below and above the normal range, respectively. Each laboratory should establish ranges with their instruments and reagents. Follow troubleshooting procedures when results fall outside the established range.

### VI. Limitations

All controls are subject to the limitations of the test system. Variables such as temperature, reagent stability, instrument performance, and individual technique can influence final results. Always follow instrument and reagent manufacturers guidelines.

### VII. Expected Values

Values were assigned using thrombin clotting time (Claus<sup>1</sup>) methodology. Recovered values may vary with other procedures. Actual values recovered depend on the instrument and reagent used.

### VIII. Performance Characteristics

Low Fibrinogen Control is formulated to have a value of approximately 70-120 mg/dL. High Fibrinogen Control should be approximately 500 to 650 mg/dL. Fibrinogen controls have been used in a quality control program covering multiple runs over a period of time. The coefficient of variation was less than 10%.<sup>4</sup>

### IX. References

1. Claus, A., Rapid physiological coagulation method for the determination of fibrinogen, *Acta Haemat* 17:237-246. 1957.
2. NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
3. Stability data found in DHF.
4. Data found in 510K File.

#### ORDERING INFORMATION

Cat. No.	Description	Contents
100604	Low Fibrinogen Control Plasma	10 x 1 mL
100607	High Fibrinogen Control Plasma	10 x 1 mL
100600	Fibrinogen Assay Set	100 determinations
100605	Fibrinogen Assay Set	300 determinations

#### FISHER DIAGNOSTICS® LIMITED WARRANTY

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

Pacific Hemostasis® is a registered trademark of Fisher Scientific Company L.L.C.  
Fisher Diagnostics® is a registered trademark of Fisher Scientific Company L.L.C.

MDCI Ltd.  
Arundel House  
1 Liverpool Gardens  
Worthing, West Sussex BN11 1SL  
UK

**Fisher Diagnostics®**  
8365 Valley Pike  
Middletown, VA 22645-0307, USA  
1-(800) 528-0494 USA only

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### IX. Bibliografia

1. Clauss, A., Rapid physiological coagulation method for the determination of fibrinogen, *Acta Haemat* 17:237-246. 1957.
2. NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and Performance of Coagulation Assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
3. I dati relativi alla stabilità sono stati ottenuti da DHF.
4. Informazioni ottenute dal file 510(k).

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