## KONTACT®

### I. Intended Use

Pacific Hemostasis® KONTACT reagent is intended for use in performing the activated partial thromboplastin time (APTT) test, and for APTT-based factor assays using a particulate activator.

### II. Summary and Principles

The APTT is useful as a screening tool, and as a quantitative test for the intrinsic coagulation factors. It is a simple and versatile test which is sensitive to deficiencies of all plasma clotting factors except Factor VII. However, it is mainly used to detect deficiencies in Factors VIII, IX, XI, XII, and Prekallikrein.

The APTT is also commonly used to monitor heparin therapy since APTT prolongation is directly proportional to increasing amounts of heparin.<sup>1,2</sup>

The APTT test is performed by adding reagent containing a plasma activator and phospholipid to the test specimen. This mixture is incubated for 5 minutes at 37°C for optimum activation. Calcium chloride is added and clot formation is timed.

### III. Reagent

For *in vitro* diagnostic use.

**Composition:** 1.2% rabbit brain phospholipid, 0.03% magnesium aluminum silica, 0.4% phenol, 0.8% buffers, salt and stabilizers

Store unopened vials at 2–8°C. **Do not freeze.** Opened vials are stable for 30 days when stored at 2–8°C. <sup>4</sup> The activator particles will settle on prolonged storage. Shake vigorously before opening to ensure adequate resuspension. Provide some mechanism such as a magnetic stirrer to maintain adequate suspension during use. Erratic values, quality control values outside established ranges, or product color variations could indicate deterioration. However, poor performance could also be due to other factors within the test system.

### IV. Specimen Collection

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22–24°C. For more details on specimen collection and storage, see NCCLS Document H21-A4.<sup>3</sup>

- Do not delay mixing the blood with anticoagulant.
- Avoid foaming the specimen.
- Use only plastic or siliconized borosilicate glass containers.
- Turbid, icteric, lipemic, or hemolyzed specimens may generate erroneous results.
- Freezing and thawing plasma that contains residual cells will generate damaged cell membranes that can affect results.
- Acute inflammatory reactions can shorten APTT results because of elevated fibrinogen.
- Plasma samples with hematocrits outside the range of 20–55% may be improperly anticoagulated and should be adjusted appropriately.

### V. Test Procedure

**Materials Provided:** KONTACT Reagent, liquid, 10 x 10 mL, or 10 x 4 mL.

**Materials Required, But Not Provided:**

Pacific Hemostasis Calcium Chloride solution (0.025M)

Stopwatch or timer

Precision pipetter: 0.1 mL

Normal and abnormal controls such as Pacific Hemostasis Coagulation Control Plasmas, Level 1, 2, and 3.

KONTACT is suitable for use with manual, mechanical, photo-optical, or other means of clot detection. Follow manufacturer's recommendations for proper use of instrumentation. For manual assays:

- Prewarm Calcium Chloride (0.025M) to 37°C.
- Add 0.1 mL test plasma to cuvette and prewarm to 37°C.
- Add 0.1 mL KONTACT to the test plasma. Mix.
- Incubate the plasma-reagent mixture at 37°C for 5 minutes (activation time).
- Forcibly add 0.1 mL prewarmed Calcium Chloride and time clot formation.

### VI. Quality Control

Normal and abnormal plasmas such as Pacific Hemostasis Coagulation Control Level 1, 2, and 3 should be tested in conjunction with patient plasmas. Level 1 is a normal plasma, and Levels 2 and 3 are adjusted to mimic moderately and severely deficient plasmas, respectively. A normal control and at least one abnormal control should be run at the initiation of testing each day and at least once each shift, or with each group of assays. Controls should also be tested with each reagent change or major instrument adjustment. Each laboratory should establish a control range to represent the allowable variation in day to day performance for each control.

### VII. Results

Report clotting times for each plasma to the nearest 0.1 second. A Normal Reference Range can also be reported for comparison. Do not report patient values relative to commercial control plasma clotting times. Controls are intended only for quality assurance of the test system.

### VIII. Limitations

The biochemistry of coagulation involves a series of reactions that are influenced by many pre-test conditions. These variables must be controlled to obtain reproducible results.<sup>3</sup>

#### Technique

- Plasma pH will increase if exposed to air. Store samples stoppered.
- KONTACT was designed to work at 37°C ± 0.5°C. Frequently check the temperature of all heating elements.
- All labware must be clean and free of trace amounts of detergents.
- Always follow instrument manufacturer's instructions for proper maintenance.

#### Interfering Substances

- Sodium oxalate, EDTA, and heparin are not suitable anticoagulants.
- Oral contraceptives, estrogen, pregnancy, coumarin type drugs, heparin, asparaginase, and naloxone have been reported to influence APTT results.<sup>4</sup>

### IX. Expected Values

When KONTACT was evaluated on a normal population, the following results were obtained: <sup>5</sup>

	mean	± 2SD range
<b>Mechanical</b>	30.3	26.1 - 34.5
<b>Photo-optical</b>	33.3	28.5 - 38.1

These values should only be used as a guideline. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent.<sup>4</sup> A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

### X. Performance Characteristics

#### Heparin Sensitivity:

The anticoagulant action of heparin depends on many factors, including an adequate level of Antithrombin-III, platelet activation and subsequent Platelet Factor-4 release during specimen preparation, *in vivo* presence of other medications, rate of heparin metabolism, mode of heparin administration, and delayed specimen handling. While recognizing these variables, the laboratory can determine the relative sensitivity of a given reagent to heparin by adding known amounts of heparin to pooled normal plasma and performing an APTT.

ENGLISH

For example, the following results were obtained on a photo-optical instrument with one lot of KONTACT reagent:<sup>10</sup>

Heparin conc. (units/mL)	APTT (seconds)
0.0	31.9
0.05	35.0
0.1	39.9
0.2	53.2
0.3	68.5
0.4	85.3
0.5	106.5

Each laboratory should establish its own heparin sensitivity curve using the same heparin source used for therapy in that institution. Variations can result from different brands of heparin, tissue origin, and salt forms.<sup>12,5</sup>

**Warning: DO NOT USE with Instrumentation Laboratory line of ACL™ COAGULATION ANALYZERS on HEPARINIZED patient samples.**

#### Factor Sensitivity:

An APTT reagent with adequate sensitivity should demonstrate a prolonged clotting time in samples having ≤ 30-40% factor activity.<sup>6,7</sup> KONTACT was evaluated on mildly and severely deficient plasmas with the following results: <sup>11</sup>

Factor	% activity	APTT (seconds)
VIII	<1%	150
VIII	20%	50.7
IX	<1%	>150
IX	20%	39.6
XI	<1%	88.7
XI	20%	44.5
XII	<1%	>150
XII	20%	43.2
Prekallikrein	<1%	55.9

Furthermore, the sensitivity of KONTACT to Factor VIII has been determined as follows: <sup>12</sup>

% Factor VIII	APTT (seconds)
100%	34.3
70%	37.8
50%	40.9
40%	43.4
30%	46.2
20%	50.7
10%	57.0
5%	63.6
1%	85.2
<1%	>150

These values should only be used as guidelines. Each laboratory should establish sensitivity to individual factors using instruments, reagents, and techniques used in their laboratory.

ORDERING INFORMATION			
	<b>Cat. No.</b>	<b>Description</b>	<b>Contents</b>
	100312	KONTACT	10 x 10 mL
	100308	KONTACT	10 x 4 mL
	100314	CaCl <sub>2</sub> (0.025M)	10 x 10 mL
	100309	CaCl <sub>2</sub> (0.025M)	10 x 4 mL

#### FISHER DIAGNOSTICS® LIMITED WARRANTY

Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD's sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

Pacific Hemostasis® is a registered trademark of Fisher Scientific Company L.L.C.

Fisher Diagnostics® is a registered trademark of Fisher Scientific Company L.L.C.

### XI. References

- Brandt, J.T, Triplett, D.A.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated Partial Thromboplastin Time. *Amer J Clin Path* 76:530, 1981.
- Thompson, J.M.: The Control of Heparin Therapy by the Activated Partial Thromboplastin Time. Sensitivity of Various Thromboplastins to Heparin. *Standardization of Coagulation Assays: An Overview*. Edited by D.A. Triplett. College of American Pathologists, Skokie, Ill, 1982, pp 195.
- NCCLS: *Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. *Clin Chem* 18:1041, 1972.
- Banez, E.L., Triplett, D.A., Koepke, J.: Laboratory Monitoring of Heparin Therapy. The Effect of Different Salts of Heparin on the Activated Partial Thromboplastin Time. *Amer J Clin Path* 74:569, 1980.
- Wujastyk, J., Triplett, D.A.: Selecting Instrumentation and Reagents for the Coagulation Laboratory. *Pathologist* 37:398, 1983.
- Christensen, R.L., Triplett, D.A.: Factor Assay (VIII and IX) Results in the College of American Pathologists Survey Program (1980-1982). *Amer J Clin Path* 80 (Suppl): 633, 1983.
- 2-12. Data found in 510(K) file.

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- Brandt, J.T, Triplett, D.A.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated Partial Thromboplastin Time. *Amer J Clin Path* 76:530, 1981.
- Thompson, J.M.: The Control of Heparin Therapy by the Activated Partial Thromboplastin Time. Sensitivity of Various Thromboplastins to Heparin. *Standardization of Coagulation Assays: An Overview*. Edited by D.A. Triplett, College of American Pathologists, Skokie, Ill. 1982, pp 195.
- NCCLS: *Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
- Young, D.S., Thomas, D.W, Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. *Clin Chem* 18:1041, 1972.
- Banez, E.L., Triplett, D.A., Koepke, J.: Laboratory Monitoring of Heparin Therapy. The Effect of Different Salts of Heparin on the Activated Partial Thromboplastin Time. *Amer J Clin Path* 74:569, 1980.
- Wujastyk, J., Triplett, D.A.: Selecting Instrumentation and Reagents for the Coagulation Laboratory. *Pathologist* 37:398, 1983.
- Christensen, R.L., Triplett, D.A.: Factor Assay (VIII and IX) Results in the College of American Pathologists Survey Program (1980-1982). *Amer J Clin Path* 80 (Suppl): 633, 1983.

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### XI. Bibliografia

- Brandt, J.T, Triplett, D.A.: Laboratory Monitoring of Heparin. Effect of Reagents and Instruments in the Activated Partial Thromboplastin Time. *Amer J Clin Path* 76:530, 1981.
- Thompson, J.M.: The Control of Heparin Therapy by the Activated Partial Thromboplastin Time. Sensitivity of Various Thromboplastins to Heparin. *Standardization of Coagulation Assays: An Overview*. Edited by D.A. Triplett, College of American Pathologists, Skokie, Ill. 1982. pp 195.
- NCCLS: *Collection, transport, and processing of blood specimens for coagulation testing and performance of coagulation assays*. 4th edition. Approved guideline. NCCLS Document H21-A4. Wayne, PA, 2003.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: Effect of Drugs on Clinical Laboratory Tests. *Clin Chem* 18:1041, 1972.
- Banez, E.L., Triplett, D.A., Koepke, J.: Laboratory Monitoring of Heparin Therapy: The Effect of Different Salts of Heparin on the Activated Partial Thromboplastin Time. *Amer J Clin Path* 74:569, 1980.
- Wujastyk, J., Triplett, D.A.: Selecting Instrumentation and Reagents for the Coagulation Laboratory. *Pathologist* 37:398, 1983.
- Christensen, R.L., Triplett, D.A.: Factor Assay (VIII and IX) Results in the College of American Pathologists Survey Program (1980-1982). *Amer J Clin Path* 80 (Suppl): 633, 1983.
- Informazioni ottenute dal file 510(k).

ITALIANO

MDCI Ltd.  
Arundel House  
1 Liverpool Gardens  
Worthing, West Sussex BN11 1SL  
Regno Unito

**Fisher Diagnostics®**  
8365 Valley Pike  
Middletown, VA 22645-0307,  
Stati Uniti  
1-(800)-528-0494 solo per gli USA

<b>Symbols Key</b>	
	Manufacturer Hersteller Fabricant Fabricante Fabbricante
	<i>In Vitro</i> Diagnostic Medical Device Medizinprodukt für die <i>in-vitro</i> -Diagnostik Matériel médical pour utilisation diagnostique <i>in vitro</i> Dispositivo médico para diagnóstico <i>in vitro</i> Dispositivo medico per diagnosi <i>in vitro</i>
	Lot Number Chargenummer Numéro de lot Número de lote Numero di lotto
	Use By Verfallsdatum Utiliser jusque Fecha de caducidad Da utilizzare entro
	Temperature Limitation Temperaturreinschränkungen Limite de température Limite de temperatura Limiti di temperatura
	CE Mark CE-Markierung Marquage CE Marca CE Marchio CE
	Catalogue Number Katalognummer Référence catalogue Número de catálogo Numero di catalogo
	Consult Instructions for Use Bedienungsanleitung lesen Consulter le manuel d'utilisation Consultar las instrucciones de uso Consultare le istruzioni per l'uso
	Pack Date Verpackungsdatum Date d'emballage Fecha del envase Data di confezionamento
	Authorized Representative in the European Community Autorisierte Vertretung in der Europäischen Gemeinschaft Représentant agréé pour La Communauté européen Representante autorizado en la Comunidad Europea Rappresentante autorizzato nella Comunità europea