Infinity™ ACE Calibrator

(Angiotensin Converting Enzyme)

INTENDED USE

The ACE Calibrator is intended for the calibration of ACE assays. It is for in vitro diagnostic use only.

PRINCIPLES OF THE PROCEDURE

This serum calibrator is designed for use as a calibrator with Infinity ACE (Angiotensin Converting Enzyme) reagent. The product consists of a lyophilized serum.

REAGENT COMPOSITION

The Calibrator contains ACE (porcine) in a buffered human serum base with stabilizers and preservatives.

WARNING: Human serum is used in the manufacture of this product. The source material used has been tested by FDA-approved methods and found non reactive for hepatitis B, surface antigen (HBsAg), antibody to hepatitis C (HCV), and the antibody to HIV1/2. These test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that hepatitis B virus (HBV), hepatitis C (HCV), and human immunodeficiency virus (HIV 1/2) or other infectious agents are absent, all products manufactured using human source material should be handled as potentially infectious. This product may also contain other human source material for which there are no approved tests. Exercise precaution in handling metal crimps and broken glass vials, as sharp edges can injure the user. For further information consult the Infinity Ace Calibrator Material Safety Data Sheet.

RECONSTITUTION

- Remove the metal seal and rubber stopper from the vial of lyophilized serum.
- Using a volumetric pipette, add exactly 1.0 mL of distilled or deionized water to the lyophilized serum.
- Replace the stopper in each vial and allow to sit for 5 to 10 minutes.
- Swirl the contents until dissolved, then remove sample for testing.

STORAGE AND STABILITY

Unreconstituted material - The ACE Calibrator should be stored between 2° to 8°C and is stable until the expiry date which appears on the kit and vial labels.

Reconstituted material - When stored capped at 2° to 8°C the constituents are stable for at least 7 days.

DETERIORATION

Excessive turbidity may indicate microbial contamination in which case the Calibrator should not be used.

PROCEDURE

The ACE Calibrator is to be used in accordance with the directions accompanying the assay procedure being used. While following these directions, the calibrator is to be treated in the same manner as a patient's specimen would be treated.

LIMITATIONS

- The results obtained using the calibrator are dependent upon several factors. Erroneous results can occur from reconstitution inaccuracies and technique errors associated with the assay procedure. For more information, refer to the "Limitations" section of the package insert of the assay procedure being used.
- Improper storage or handling of the calibrator can also affect results. If there is visible evidence of microbial growth in a vial, do not use that vial or product.

ASSAY VALUE

Please refer to vial label, for calibrator value.

The assigned value for the Infinity ACE Calibrator, traceable to an in-house master calibrator, has been established using current lots of Infinity ACE reagents. For use with other brands of reagents the user should assure themselves that the target values are still valid. The procedure of Miller WG, et al., Clin. Chem 1982; 28:2195-2200 is recommended for this purpose.

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© 2009 Thermo Fisher Scientific Inc. All rights reserved.