Infinity™ ***Ammonia Standard Set***

INTENDED USE

This standard set is intended for the calibration and control of Ammonia assays using the Thermo Infinity Ammonia reagent (TR60101). It is for in vitro diagnostic use only.

PRODUCT COMPOSITION

The Thermo Infinity Ammonia standard set is prepared from a stock Ammonium salt solution, with ammonia free water, and are supplied as follows:

LEVEL 1	59 µmol/L (100 µg/dL) 1 mg/L Ammonium Salt
LEVEL 2	118 µmol/L (200 µg/dL) 2 mg/L Ammonium Salt
LEVEL 3	177 µmol/L (300 µg/dL) 3 mg/L Ammonium Salt

WARNING: Do not ingest. Avoid contact with skin and eyes. If spilt thoroughly wash affected areas with water. Standard contains sodium azide which may react with copper or lead plumbing. Flush with plenty of water when disposing. For further information consult the Thermo Infinity Ammonia Standard Set Material Safety Data Sheet.

PREPARATION

The standards are provided ready to use.

STORAGE AND STABILITY

The standards are stable until the expiry date stated on the bottle label when stored at 10-25°C. After opening, the bottles should be kept tightly capped to prevent evaporation.

The solutions should be clear, colourless and free of precipitate. Any discolouration or turbidity indicates deterioration and the solution should be discarded.

PROCEDURE

The Thermo Infinity Ammonia Standard Set should be handled in the same manner as a patient sample following the assay procedure detailed in the Thermo Ammonia Reagent package insert.

PERFORMANCE

Each batch of standards is assayed using Thermo Infinity Ammonia reagent to confirm the stated value.

LIMITATIONS

These standards have been validated for use with the Thermo Infinity Ammonia Reagent (TR60101). Users should verify its suitability for use with other Ammonia methods. Always dispense standard using the dispenser nozzle into a clean analyser cup. Do not place pipettes into the standard bottle. Keep capped when not in use.

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© 2008 Thermo Fisher Scientific Inc. All rights reserved.
a division of Fisher Scientific Company, LLC
a subsidiary of Thermo Fisher Scientific Inc.
MDCI Ltd.