Chloride Standard

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Chlorid-Standard

INTENDED USE

This standard is intended for the calibration of assays for the quantitative determination of Chloride. It is for in vitro diagnostic use only.

REAGENT COMPOSITION

The standard is an aqueous solution containing 100 mEq/L (100 mmol/L) of Chloride, with an added preservative.

WARNING

Do not ingest. Avoid contact with skin and eyes. If spilt, thoroughly wash the affected area with water.

PREPARATION

The standard is ready to use as supplied.

STORAGE AND STABILITY

The unopened bottles are stable to the expiry date stated on the label when stored at 2-25 °C. After opening, the bottles should be kept tightly capped to prevent evaporation. The solution should be clear and colourless. Any discolouration or turbidity indicates deterioration and the solution should be discarded.

PROCEDURE

The Chloride Standard is suitable for use as a calibrator with most commonly used colourimetric methods for Chloride determination. The standard should be handled in the same manner as a patient sample following the assay procedure detailed in the package insert of the reagent being used.

PERFORMANCE

Each batch of standard is assayed using Thermo colourimetric Chloride reagent to confirm the stated value.

TRACEABILITY

The standard has been verified against available Standard Reference Materials from the US National Institute of Standards and Technology (NIST).

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© 2008 Thermo Fisher Scientific Inc. All rights reserved.
a division of Fisher Scientific Company, LLC
a subsidiary of Thermo Fisher Scientific Inc.
MDCI Ltd.