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## Oxalate Urine Controls Normal and Elevated

### INTENDED USE

The Oxalate Normal and Elevated Controls are for use in the quality control of the Thermo Scientific Oxalate method for the quantitative determination of urine oxalate. The control is used for monitoring accuracy and precision on automated clinical chemistry analysers and for manual operations. It is for in vitro diagnostic use only.

### PRODUCT COMPOSITION

The Oxalate Normal and Elevated Controls are aqueous synthetic urine preparations, formulated for monitoring the accuracy and precision of urine oxalate assays.

**WARNING:** Do not ingest. Avoid contact with skin and eyes. If spilt thoroughly wash affected areas with water. For further information consult the Oxalate Normal and Elevated Controls Material Safety Data Sheet.

### PREPARATION

The controls are provided ready to use.

### STORAGE AND STABILITY

The controls are stable until the expiry date stated on the bottle label when stored at 2-8°C. After opening, the bottle should be kept tightly capped to prevent evaporation.

### DETERIORATION

The solution should be clear, faint yellow and free of precipitate. Any turbidity indicates deterioration and the solution should be discarded.

### PROCEDURE

The Thermo Scientific Oxalate Normal and Elevated Controls should be handled in the same manner as a patient sample following the assay procedure detailed in the Thermo Scientific Oxalate Reagent package insert.

### LIMITATIONS

- A variety of factors may affect the results obtained using the Oxalate Normal and Elevated Controls. Among these factors are improper pipetting, inadequate mixing and technique errors associated with the particular instrument or reagent system.
- The Oxalate Normal and Elevated Controls should not be used past the expiration date.

### TARGET VALUE AND ACCEPTABLE RANGES

The following table provides target values and acceptable ranges for Oxalate Urine Controls.

<b>Control</b>	<b>Target Value</b>	<b>Acceptable Range</b>
Normal	0.30 mmol/L	0.25 – 0.35 mmol/L
	27 mg/L	22.5 – 31.5 mg/L
Elevated	1.00 mmol/L	0.84 – 1.16 mmol/L
	90 mg/L	75.6 – 104.4 mg/L

The target values were verified using Thermo Scientific Oxalate reagents. Determinations were performed using current lots of Thermo Scientific reagents on the Roche Cobas Mira, using approved standardised in-house procedures. The target values are traceable to NIST reference material. The assigned or acceptable range, is the expected range of results to be used as a guide in determining the accuracy of the procedure in each laboratory. Because of differences existing between instrumentation, techniques and reagents, the mean value achieved by a laboratory should fall within the corresponding acceptable range. It is recommended that each laboratory establish its own mean and acceptable range and use those provided only as a guide.

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© 2008 Thermo Fisher Scientific Inc. All rights reserved.