Oxalate Standardprøve

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Oxalate Standard

INTENDED USE

This standard is intended for the calibration of Oxalate assays using the Thermo Scientific Oxalate Reagent (TR86001). It is for in vitro diagnostic use only.

PRODUCT COMPOSITION

The Thermo Scientific Oxalate Standard is prepared gravimetrically from a certified primary reference material, and is supplied as an aqueous preparation with preservative. The standard contains 0.50 mmol/L Oxalate.

WARNING: Do not ingest. Avoid contact with skin and eyes. If spilt thoroughly wash affected areas with water. For further information consult the Thermo Scientific Oxalate Standard Material Safety Data Sheet.

PREPARATION

The standard is provided ready to use.

STORAGE AND STABILITY

The standard is stable until the expiry date stated on the bottle label when stored at 2-8°C After opening, the bottle should be kept tightly capped to prevent evaporation.

DETERIORATION

The solution should be clear, colourless and free of precipitate. Any discolouration or turbidity indicates deterioration and the solution should be discarded.

PROCEDURE

The Thermo Scientific Oxalate Standard should be used neat, without the sample pre-treatment protocol provided for Thermo Scientific Oxalate Controls and patient samples. For further details, refer to the assay procedure detailed in the Thermo Scientific Oxalate Reagent package insert.

PERFORMANCE

Each batch of standard is assayed using Thermo Scientific Oxalate reagent to confirm the stated value.

LIMITATIONS

This standard has been validated for use with the Thermo Scientific Oxalate Reagent (TR86001). Users should verify its suitability for use with other Oxalate methods. Always dispense standard using the dispenser nozzle into a clean analyser cup. Do not place pipettes into the standard bottle. Keep capped when not in use.

TRACEABILITY

Thermo Scientific Oxalate Standard is traceable to the US National Institute of Standards and Technology (NIST) reference material.

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© 2008 Thermo Fisher Scientific Inc. All rights reserved.