












PBG Test Kit

SYMBOLS IN PRODUCT LABELLING

	Authorised Representative		Temperature Limitation
	For in vitro diagnostic use		Use by/Expiration Date
	Batch code/Lot number		CAUTION. Consult instructions for use.
	Catalogue number		Manufactured by
	Consult instructions for use		Xn - Harmful
			C - Corrosive

INTENDED USE

This kit is intended for the in vitro qualitative determination of porphobilinogen (PBG) in urine.

CLINICAL SIGNIFICANCE^{1,2,3,4}

The porphyrias are a group of disorders which result from abnormalities in the biosynthesis of haem. They can be divided into two groups, the acute porphyrias and the non-acute (predominantly cutaneous) porphyrias. A number of clinically distinct disorders occur in each group. The three common acute porphyrias are:

- (i) Acute Intermittent Porphyria (AIP)
- (ii) Porphyria Variegata (PV); and
- (iii) Hereditary Coproporphyrinuria (HC).

AIP, PV and HC are dominantly inherited disorders. They may present in an acute phase (with neurological symptoms) or in a latent phase. PV and HC may have skin symptoms in addition to neurological symptoms.

The presenting features of an acute attack include abdominal pain and neurological symptoms ranging from peripheral neuritis to quadriplegia. If an acute attack is not diagnosed, the patient may be subjected to surgery, with the use of anaesthetics, which may further aggravate the condition. In severe attacks weakness of the trunk muscles can cause respiratory failure and sometimes death.

Several factors, including exposure to a variety of common drugs, changes in hormonal status, diet or acute illness may precipitate an acute attack.

During an acute attack the haem precursor, porphobilinogen (PBG) accumulates in the liver, and raised levels occur in the plasma and urine. An increase in urinary PBG therefore is strongly indicative of an acute porphyria attack.

METHODOLOGY^{4,5,6}

The PBG screening test kit is based on the Watson - Schwartz test. In the Watson-Schwartz test, PBG present in the urine condenses with p-dimethylaminobenzaldehyde (DMAB) in acid solution to form a magenta coloured product. In order to improve the specificity of the method it is necessary to then remove substances which may interfere. The most common one being urobilinogen which produces a magenta colour similar to PBG with DMAB/acid solution. Organic extractions with centrifugations are usually carried out to remove these interfering substances. These extractions are time consuming and a certain level of expertise is required to interpret the test.

The PBG screening method utilises an anion exchange resin which binds PBG present in the urine. Interfering compounds are then removed through a simple washing step. PBG is then eluted off the resin and added to the DMAB/acid solution. If PBG is present in the sample in abnormal amounts a magenta colour develops. By comparing the colour developed with the colour chart included in the kit, an approximate concentration of PBG can be determined.

Positive results should be confirmed by quantitative estimation. Quantitative tests require a high level of expertise and are best carried out by a reference laboratory specialising in the investigations of porphyria.

KIT COMPOSITION

ACTIVE INGREDIENTS	QUANTITY
Resin filled syringes	20
Filters (5 micron)	20
Dimethyl aminobenzaldehyde Powder(DMAB)	2x10 Tests (10mL)
DMAB Diluent (6.1 Molar HCl)	2x10 Tests (10mL)
Reaction tubes	20
Elution Reagent (1 Molar Acetic Acid)	20 Tests (20mL)
Colour Chart	1

WARNING:

- R22 Harmful if swallowed.
R34 Causes burns.
R36/38 Irritating to eyes and skin.
R37 Irritating to respiratory system.
S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S28 After contact with skin, wash immediately with plenty of soap and water.

For further information, consult the PBG test kit Material Safety Data Sheet.

REAGENT PREPARATION

Reconstitute the DMAB with DMAB diluent as follows:-

1. Remove cap from DMAB diluent and place pouring adapter onto the neck of the bottle.
2. Remove cap from DMAB and pour the contents of the DMAB diluent into the DMAB bottle. Mix gently.
3. Discard the Diluent bottle, cap and adapter.
4. When stored at room temperature (18-25°C) in the bottle supplied, this reagent is stable for at least 12 months. Record the expiry date (reconstitution date plus 12 months) in the place provided on the DMAB label.

STABILITY AND STORAGE

All reagents are stable until the expiry date shown on the label when stored at room temperature (18-25°C).

SPECIMEN COLLECTION AND HANDLING⁴

Urine: A random urine sample is suitable for use with this method. Urine should be collected without the use of a preservative or stabiliser.

Storage: The urine sample should be protected from light at all times. Urine should be tested for PBG within 8 hours of voiding. If this is not possible then store at -20°C or less. Samples stored at -20°C and protected from the light are stable for at least 6 months.

ADDITIONAL EQUIPMENT REQUIRED

- pH meter, pH paper or pH indicator solutions.
- Pipette.
- Distilled or deionised water.
- Timer.
- Ammonia solution, approximately 8% (Concentrated ammonia solution is approx. 33%. Dilute this 1:4 with distilled water to approx. 8%).

