Shandon MUCOLEXX™ Mucoliquefying Preservative

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

Shandon MUCOLEXX™

Preservation and Mucoliquefying Agent

MUCOLEXX is a mucoliquefying preservative designed for use in the preparation of mucoid cytology specimens. For years cytotechnologists have been plagued with the problem of sticky, thick mucoid specimens. Mucolexx liquefies mucus and allows for optimal cytologic diagnostic results.

WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE.

Wash hands thoroughly after use.
Keep container closed.

STORAGE

Keep tightly closed. The presence of a small amount of bottom precipitate does not indicate a breakdown of the solution. Store at room temperature, less than 90° F (32° C). No refrigeration is required. Protect from exposure to chemical fumes. Shelf life is two years unopened. After opening, keep container tightly sealed and stored properly. To ensure reagent effectiveness, mix fresh cells (e.g., blood cells) with a small amount of MUCOLEXX on a clean glass slide. When examined microscopically, the cells should not appear distorted due to shrinkage or swelling (i.e. hypo- or hypertonicity).

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

MIXING INSTRUCTIONS

Use full strength – do not dilute.

SPECIMEN PREPARATION

The specimen may be collected in a cup which contains 5-10 times the specimen volume of MUCOLEXX™. Pre-filled disposable cups are available from Thermo Electron Anatomical Pathology. This will preserve the specimen until it reaches the laboratory. In cases where this would be inconvenient, use a clean, dry collection cup. **DO NOT** collect in alcohol. Specimen should be vigorously shaken as soon as possible after collection in MUCOLEXX. Use the following method for fresh specimens:

1. Place fresh, unfixed specimen in a 50 ml test tube.
2. Add an equal volume of MUCOLEXX to the specimen. Depending upon the degree of specimen viscosity, additional MUCOLEXX may be required. **CAUTION: MUCOLEXX IS NOT A STERILIZING AGENT.** When handling specimens in MUCOLEXX, universal precautions should be followed to avoid possible infection.
3. Cap the test tube and shake vigorously by hand. This step usually takes 30-60 seconds.
4. At this point the specimen may be centrifuged, the supernatant decanted, and the resulting button used to prepare smears, Cytospin® preparations, and/or cell blocks.

RESULTS

Mucus is liquefied. Microscopically, cells exhibit no cellular distortion or destruction. Cells display normal staining characteristics.

LIMITATIONS

It has been found that some specimens do not liquefy readily. If additional MUCOLEXX, additional time and additional agitation do not liquefy the mucus, one can disperse mucus by use of a blender. Cells do not tend to adhere to the blender when MUCOLEXX is used. For certain tenacious mucus samples which do not liquefy adequately in MUCOLEXX, the Saccomanno Blender technique may be preferred. Pre-filled collection cups of fluid suitable for this use are available from Thermo Electron Anatomical Pathology.

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