# Shandon Rapid-Chrome™ Papanicolaou Staining Kit

**Thermo**  
ELECTRON CORPORATION

Anatomical Pathology USA  
Clinical Diagnostics  
171 Industry Drive  
Pittsburgh, PA 15275, USA  
Tel: 1-800-547-7429  
+1 412 788 1133  
Fax: +1 412 788 1138  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

Anatomical Pathology International  
Clinical Diagnostics  
93-96 Chadwick Road  
Astmoor, Runcorn  
Cheshire, WA7 1PR, UK  
Tel: +44 (0) 1928 566611  
Fax: +44 (0) 1928 565845  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

## SHANDON RAPID-CHROME™ PAPANICOLAOU STAINING KIT

This staining kit is specially designed to facilitate staining of cytology specimens. Fresh specimens can be stained within minutes of their collection.

With this Rapid-Chrome™ Staining Kit, all the reagents for permanently stained cytology specimens are contained in a portable tray. The standard Papanicolaou stains are used: Gill 2 Hematoxylin, EA-50 and Orange G-6. All reagents are sequentially aligned to allow the slide to be moved from one jar to the next. Up to 100 slides can be stained with this kit.

The Rapid-Chrome™ Papanicolaou Staining Kit is ideally suited to quickly stain fine needle aspirates. It also prevents cross contamination of routine staining dishes from high-risk cases.

### WARNINGS AND PRECAUTIONS

#### FOR IN VITRO DIAGNOSTIC USE.

#### **WARNING: CONTAINS FLAMMABLE LIQUIDS.**

Reasonable care should be taken when using all laboratory reagents. Use only with adequate ventilation. Keep away from direct sunlight.

---

### WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

---

**NOTE:** This kit contains Reagent Alcohol (100% and 95%), Xylene, and Rinse water, which are considered general laboratory reagents.

### STORAGE

Store at room temperature. Do not handle or store near heat, sparks, flames, or strong oxidants.

Lot number and expiration date are printed on the tray.

## INSTRUCTIONS FOR USE

### A. Kit Preparation

1. All solutions required are provided in the tray and are ready to use full strength.
2. Use distilled water in the Hematoxylin jar to replace water lost by evaporation.
3. Bluing Reagent should be changed as often as necessary to retain its alkalinity. A 500 ml replacement bottle is supplied with each kit for this purpose.
4. Other solutions requiring periodic replacement are:
  - Distilled Water – frequently
  - Alcohols – as needed (ethanol or reagent)
  - Xylene – as needed

### B. Specimen Preparation

1. Cytology specimens must not be allowed to air dry prior to fixation. They should be fixed immediately in 95% alcohol or with a spray fixative such as Shandon Cell-Fixx™, or in a collection fluid fixative such as Shandon Cytospin® Collection Fluid.
2. Fresh cytology specimen may be placed directly into the first 95% alcohol jar in the Rapid-Chrome Staining Kit.
3. Specimens which were prepared with a carbowax fixative require adequate rinsing to remove the coating before staining. Failure to completely remove the coating will result in understaining. If you use a commercial fixative, follow the manufacturer's recommendation for removal.

### C. Staining Procedure

1. Remove all lids and place each in its slot to the right of the jar.
2. Place slide in the jar of 95% alcohol next to the "start" arrow (lower left corner of tray).
3. Process slides through the solutions following the directions indicated by the arrows.
4. A recommended staining protocol can be found on the following page. Staining times may be adjusted to suit individual preference in stain intensity.
5. Coverslip the slide using Shandon Mounting Medium.

**NOTE:** Thermo has not validated the protocol illustrated here and takes no responsibility for its use. Customers should always validate protocols before placing any reliance on them.

## RECOMMENDED RAPID-CHROME™ PAPANICOLAOU STAIN PROTOCOL

Step	Solution	Time
1	95% Alcohol	1 Minute
2	95% Alcohol	10 Dips
3	Distilled Water	10 Dips
4	Hematoxylin Gill 2	1 Minute
5	Distilled Water	10 Dips
6	Bluing Reagent	1 Minute
7	Distilled Water	10 Dips
8	95% Alcohol	10 Dips
9	Orange-G-6	1 Minute
10	95% Alcohol	10 Dips
11	95% Alcohol	10 Dips
12	EA-50	1 Minute
13	95% Alcohol	10 Dips
14	95% Alcohol	10 Dips
15	100% Alcohol	10 Dips
16	100% Alcohol	10 Dips
17	Xylene	10 Dips
18	Xylene	1 Minute

**ACTIVE INGREDIENTS**

## 100% Alcohol:

Specially Denatured Alcohol Formula 3.A	95 Parts by volume
Isopropyl Alcohol	5 Parts by volume

DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.  
CANNOT BE MADE NONPOISONOUS. FLAMMABLE.

## 95% Alcohol:

Reagent Alcohol

DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.  
CANNOT BE MADE NONPOISONOUS. FLAMMABLE.

## Gill 2 Hematoxylin:

Aluminum sulfate	3.3%
Hematoxylin	0.376%

WARNING: HARMFUL OR FATAL IF SWALLOWED.

## Bluing Reagent:

Sodium bicarbonate	0.18%
Lithium carbonate	0.85%

DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.  
CANNOT BE MADE NONPOISONOUS. FLAMMABLE.

## Orange-G-6:

Orange-G	0.67%
Phosphotungstic acid	0.02%

DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.  
CANNOT BE MADE NONPOISONOUS. FLAMMABLE.

## EA-50:

Eosin-Y	0.45%
Phosphotungstic acid	0.2%
Lt. Green SF Yellowish	0.1%

DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED.  
CANNOT BE MADE NONPOISONOUS. FLAMMABLE.

## Xylene:

Xylene	100%
--------	------

DANGER. HARMFUL OR FATAL IF SWALLOWED. IF SWALLOWED, DO NOT INDUCE VOMITTING.  
CALL A PHYSICIAN IMMEDIATELY. VAPOR HARMFUL. FLAMMABLE.

## Shandon Mounting Medium:

Methacrylate polymer	30.8%
----------------------	-------

DANGER. HARMFUL OR FATAL IF SWALLOWED. IF SWALLOWED, DO NOT INDUCE VOMITTING.  
CALL A PHYSICIAN IMMEDIATELY. VAPOR HARMFUL. FLAMMABLE.

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