Shandon M-1 Embedding Matrix

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

Shandon M-1 Embedding Matrix™ is an outstanding embedding medium used to support tissue for frozen sectioning at temperatures between -12° to -14° C. It forms a firm bond with the specimen holder that secures tissue during microtomy. M-1 Embedding Matrix freezes quickly, sections smoothly, and dissolves in water without residue.

SHELF LIFE

Shandon M-1 Embedding Matrix has an unopened shelf life of 2 years from date of manufacture. Once opened, it should remain usable for 6 months if kept capped and properly stored.

STORAGE

Shandon M-1 Embedding Matrix should be stored at room temperature. Keep containers tightly closed when not in use.

FOR IN VITRO DIAGNOSTIC USE.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions.

Definitions are at the rear of this booklet

INSTRUCTIONS FOR USE:

1. Cover the surface of the specimen holder with a thick layer of M-1 Embedding Matrix™.
2. Position and orient tissue on holder, and press lightly to eliminate air bubbles beneath the tissue.
3. Allow tissue and Matrix to slightly begin freezing and then apply additional M-1 Embedding Matrix until the tissue is completely surrounded.
4. Freeze tissue to holder in usual manner until tissue and M-1 Embedding Matrix becomes firm.
5. Section specimen in usual manner.
6. No special procedure is required to remove M-1 Embedding Matrix from sections; it will dissolve during staining.
7. To remove specimen from holder for further processing, dissolve M-1 Embedding Matrix in water and immediately place specimen in fixative.

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