Shandon Staining Products

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

Shandon Staining Products provide a total solution to your staining needs.

Shandon stains are intended for use in histological and cytological staining procedures. All are formulated and packaged to provide our customers with the best value. Every reagent and raw material is tested to meet the highest standards. Only those dyes certified by the Biological Stain Commission are used in manufacturing. Each lot of stain is subjected to a series of quantitative and qualitative tests. These tests ensure maximum performance, stability and lot to lot consistency. From scientific testing and spectrophotometry, to staining a simple control sample, we make every effort to guarantee your satisfaction with Shandon stains and reagents.

Shandon Gill Hematoxylins offer speed, convenience and the choice of stain intensity. Gill 1 is recommended for cytology preparations. For darker nuclei or shorter staining times Gill 2 may be used. Gill 2 or Gill 3 is recommended for routine histology, immunocytochemistry and frozen section applications. Gill 3 is also recommended for the staining of resin embedded tissue sections.

Harris Hematoxylins are ready-to-use in your choice of acidified or non-acidified formulations. The Shandon brand does not employ a mercury-based oxidizing agent, thus removing a possible source of environmental contamination.

Shandon Bluing Reagent is a gentle, pH-controlled solution for bluing hematoxylin. It replaces the strong alkaline rinses or long tap water washes which may cause loss of cells or tissue sections.

Nu-Clear™ improves cell clarity by removing excess hematoxylin from the tissue in regressive staining protocols. It may also be used to remove background hematoxylin from the slide. Nu-Clear is designed for hand staining as well as use on automatic stainers. It is available in 2 concentrations: Nu-Clear I is faster acting, with a 10-second differentiating step; Nu-Clear II will accomplish similar results in 20 seconds. The less concentrated Nu-Clear II is recommended for use on Linear stainers and with procedures that utilize Gill Hematoxylin. The usage and timing depend on the staining intensity preferred by the laboratory.

Cytoplasmic Counterstains are offered in a complete range for use with routine H & E or Papanicolaou techniques. The range includes Instant Eosin and alcoholic or aqueous Eosin Y for Histology, as well as EA-50, EA-65 and OG-6 for Cytology. All stains produce crisp, brilliant staining to compliment the nucleus and differentiate various cellular components.

Shandon Xylene Substitute is an odorless Xylene Substitute. A safe and convenient alternative to xylene it is an aliphatic hydrocarbon that is not greasy. It may be used on stainers, automatic coverslipppers and tissue processors.

GENERAL COMMENTS

1. Shandon stains are designed to be used full strength.
2. If less intense cytoplasmic staining is desired, Eosin-Y can be diluted with 80-85% alcohol (ethanol reagent or alcohol).
3. Harris Hematoxylin and Gill 3 should be filtered before use.

4. The coating agent found in Cell-Fixx™ and Cytospin® Collection Fluid must be removed before staining. Generally, soaking in 95% alcohol will remove the coating. Failure to completely remove the coating will result in poor and uneven staining. Users of commercial spray fixatives should follow the manufacturer's recommendation for removal.
5. If using a mercuric chloride fixative (i.e. B-5, Zenkers) for histology specimens remove the mercuric chloride crystals with iodine and clear with sodium thiosulfate after removing the wax and before the hematoxylin step.
6. The pH of tap water varies considerably. Use of tap water in the staining procedures may give unpredictable results. Long-running tap water rinses are not recommended. Chlorine in tap water will bleach out the stain.
7. Guidelines for cytology and histology staining procedures can be found at the rear of the booklet. Suggested staining times and number of rinses can be adjusted to suit individual preferences in stain intensity.

NOTE: Thermo has not validated the protocols illustrated on page ____ and takes no responsibility for their use. Customers should always validate protocols before placing any reliance on them

WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE. Reasonable care should be taken when using all laboratory reagents. Use only with adequate ventilation. Keep away from direct sunlight. Contains FLAMMABLE liquids.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet.

STORAGE

Store away from direct sunlight and at room temperature. Do not freeze. Do not handle or store near heat, sparks, flames, or strong oxidants.

Histology Staining Guidelines			Cytology Staining Guidelines		
STEP	SOLUTION	TIME	STEP	SOLUTION	TIME
1	Shandon Xylene Substitute or Xylene	3-4 minutes	1	95% Alcohol (if needed)	10-15 minute pre-soak
2	Shandon Xylene Substitute or Xylene	3-4 minutes		Distilled Water Wash	1 minute
3	Shandon Xylene Substitute or Xylene	3-4 minutes	2	Harris Hematoxylin (filtered immediately before use)	.3 minutes
4	Alcohol, 100%	1 minute		OR Gill Hematoxylin (depending on formulation 1, 2 or 3)	1 to 4 minutes
5	Alcohol, 95%	1 minute	3	Distilled Water Wash	30-60 seconds
6	Alcohol (optional)	1 minute	4	Nu-Clear I	10 seconds
7	Distilled Water Wash	1-2 minutes		OR Nu-Clear II	20 seconds
8	Harris Hematoxylin (filter before use)	4-8 minutes	5	Distilled Water Wash	30-60 seconds
	OR Gill Hematoxylin (depending on formulation 2 or 3)	2-5 minutes	6	Bluing Reagent	1 minute
9	Distilled Water Wash	1-2 minutes	7	Distilled Water Wash	30-60 seconds
10	Nu-Clear I	10 seconds	8	95% Alcohol	10 seconds
	OR Nu-Clear II	20 seconds	9	Orange G-6	1-3 minutes
11	Distilled Water Wash	1 minute	10	95% Alcohol	10 seconds
12	Bluing Reagent	1 minute	11	95% Alcohol	10 seconds
13	Distilled Water Wash	1-2 minutes	12	EA-50	3 to 6 minutes
14	Alcohol (70%-95% optional)	1 minute		OR EA-65	3 to 6 minutes
15	Eosin Y, alcoholic	10 seconds to 1 minute	13	95% Alcohol	10 seconds
	OR Eosin, aqueous	2-4 minutes	14	95% Alcohol	10 seconds
16	Alcohol, 95%	20-30 seconds	15	Absolute Alcohol	20 seconds
17	Absolute Alcohol	1-2 minutes	16	Absolute Alcohol	20 seconds
18	Absolute Alcohol	1-2 minutes	17	Shandon Xylene Substitute or Xylene	1 minute
19	Absolute Alcohol	1-2 minutes	18	Shandon Xylene Substitute or Xylene	1 minute
20	Shandon Xylene Substitute or Xylene	1-2 minutes	19	Shandon Xylene Substitute or Xylene	5 minutes
21	Shandon Xylene Substitute or Xylene	1-2 minutes	20	Mount with Histo-Mount or Shandon-Mount	
22	Shandon Xylene Substitute or Xylene	1-2 minutes			
23	Mount with Histo-Mount or Shandon-Mount				

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Use by

Protect from
Heat

Thermo Electron Corporation Anatomical Pathology has realized since its beginnings the value of in-house assistance. If you have questions about ordering, shipping, delivery, inventory status or pricing, just phone our Customer Service Department toll-free at 1-800-547-7429 in the U.S., or our direct subsidiaries in the UK, France, Germany and Australia. Our technical support staff consists of certified Cytotechnologists and Histologists.