Shandon Glyo-Fixx™ Ready-to-use
Shandon Glyo-Fixx™ Concentrate

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

GLYO-FIXX™

GLYO-FIXX is formulated for use as a histological fixative, to replace formalin and zinc formalin for routine specimens as well as those intended for immunohistochemical or special staining procedures. It is also effective as a secondary fixative after unbuffered or buffered formalin or zinc formalin. GLYO-FIXX is entirely compatible with other common fixatives; phosphate and zinc salts are stable in GLYO-FIXX and will not precipitate.

The active ingredient in GLYO-FIXX is glyoxal, a two-carbon dialdehyde that fixes tissues by binding to macromolecules having methyl, amide, amine and hydroxyl groups. While some cross-linking may occur, the linkages are longer than with formaldehyde and immunoreactivity of tissue antigens is not impaired. Therefore epitope enhancement procedures such as enzyme digestion and heat treatment are often not necessary for tissues fixed in GLYO-FIXX (Initial fixation in formalin without zinc is not recommended for samples intended for Immunohistochemistry; secondary fixation in GLYO-FIXX will not reverse the effects of formalin).

There is not enough alcohol in GLYO-FIXX to have an alcoholic fixative effect. The alcohol content in ready-to-use GLYO-FIXX serves only to increase the effect of glyoxal. If you need an alcoholic fixative, GLYO-FIXX Concentrate may be diluted with alcohol (see Instructions for Use).

WARNINGS AND PRECAUTIONS FOR IN VITRO DIAGNOSTIC USE.

WARNING: FLAMMABLE

Storage

Store at room temperature. Keep container tightly closed when not in use. Store in a flammable liquid storage area.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

INSTRUCTIONS FOR USE

1. GLYO-FIXX™ Ready-to-Use is to be used full strength.
2. GLYO-FIXX Concentrate needs to be diluted.
 - a. For regular (non-alcoholic) GLYO-FIXX, mix:
 - 1 gallon (3.8 liters) GLYO-FIXX concentrate
 - 11.6 liters of deionized or distilled water
 - 3.5 liters of (100%) ethanol or reagent alcohol
 - WARNING: DO NOT USE ISOPROPANOL**
 - b. To make alcoholic GLYO-FIXX, mix:
 - 1 gallon (3.8 liters) GLYO-FIXX Concentrate
 - 7.7 liters of deionized or distilled water
 - 7.4 liters of (100%) ethanol or reagent alcohol
 - WARNING: DO NOT USE ISOPROPANOL**
3. Small specimens, such as needle biopsies, should be fixed for at least 4 hours in GLYO-FIXX. General surgical specimens require at least 6 hours after grossing, and larger specimens require longer periods.
4. For optimal results, use GLYO-FIXX for initial tissue fixation, for holding specimens after grossing, and in the first two stations on the tissue processor. To help process fatty specimens such as colon and breast, use alcoholic GLYO-FIXX in station 2 of the tissue processor.
5. Specimens may remain in GLYO-FIXX for several weeks without detriment to staining or immunoreactivity.

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