Shandon Immu-Mount™
Aqueous Non-fluorescing Mounting Medium

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

Shandon Immu-Mount™ is a water-based mounting medium specially formulated for use with frozen section, immunoperoxidase techniques, and immunofluorescent preparations.

AEC and chloro-naphtol procedures yield end-products, which are soluble in alcohol and organic solvents. Slides which have been subjected to these methods should be coverslipped with an aqueous mountant such as Immu-Mount. With proper technique, there will be no stain fading for at least one year when using Immu-Mount.

Immunofluorescent techniques require a mounting medium which will not interfere with test results. When used with filters appropriate for fluorescein and rhodamine conjugates, Immu-Mount will not exhibit background fluorescence. In addition, there is no danger of quenching fluorescent emission.

GENERAL COMMENTS

FOR IN VITRO DIAGNOSTIC USE.

Reasonable care should be taken when using all laboratory reagents. Use with adequate ventilation. Store at room temperature. Keep container closed when not in use.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

INSTRUCTIONS FOR USE

1. Remove cap and pierce hole in dispensing nozzle.
2. After saline rinse, shake excess liquid from slide and wipe back of slide.
3. Invert bottle to displace bubbles from nozzle. This will ensure a continuous column of liquid from the nozzle. **DO NOT** shake the bottle.
4. Squeeze out the first drop onto a paper towel to assure removal of air from the nozzle.
5. Without releasing pressure from the sides of the bottle, deposit several drops of mounting medium to the slide.
6. Carefully place coverslip on slide.
7. Immu-Mount may be applied to the coverslip rather than the slide: Move the nozzle along the edge of the coverslip to deposit a continuous bead of mounting medium. Place the slide at a 90° angle to the coverslip and lower the slide to a lesser angle until the coverslip adheres to the slide.
8. Use forceps to tease out any corner bubbles caught under the coverslip. If this cannot be done, soak the slide in saline, remove the coverslip, and repeat the coverslipping procedure.
9. The slide may be viewed immediately, taking care not to dislodge the coverslip before the slide has dried.
10. Allow sufficient time for the slides to dry before storage.
11. For maximum retention of fluorescence, finished slides should be stored in a cold (4° C), dark place.

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