Shandon Zinc Formal-Fixx™
Ready-to-use and Concentrate
10% Zinc Formalin
Phosphate-Free
Mercury-Free

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

ZINC FORMAL-FIXX™

Zinc Formal-Fixx Ready-to-Use and Zinc Formal-Fixx Concentrate are designed to be used in the laboratory as a routine fixative, denaturing tissue specimens and achieving cellular rigidity without over-hardening. Zinc Formal-Fixx Ready-to-Use retards protein crosslinking responsible for masking the immuno-chemistry antigenic binding sites. This fixative will yield excellent results in routine H & E, special stains enzymatic stains and immunochemical reactions.

WARNINGS AND PRECAUTIONS FOR IN VITRO DIAGNOSTIC USE

Reasonable care should be taken when using all laboratory reagents. Use formalin fixatives with adequate ventilation.

Storage:

Store at room temperature. Keep containers tightly closed. **DO NOT** store with hydrochloric acid.

Disposal:

Dispose in accordance to state and local codes. Disposal may be easier if zinc is removed from solution. Simply raise the pH to 7.0 – 8.0 with a mild alkaline solution (sodium hydroxide is recommended). Wait several hours and filter the solution. Dispose of the filtered precipitate as you would any inorganic solids of low toxicity.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

INSTRUCTIONS FOR USE

1. Zinc Formal-Fixx Ready-to-Use is to be used full strength.
2. Zinc Formal-Fixx Concentrate needs to be diluted. Mix 1 part of the concentrate with 4 parts distilled or deionized water. A one liter bottle of concentrate will make 5 liters of 10% neutral buffered formalin.
3. Use as routine histologic fixatives. For optimum results, the tissue-to-fixative ratio should range from at least 1:15 to 1:20. Specimens should be placed into fixative immediately following surgical resection for a period of time that will ensure adequate fixation.
4. It may also be used as the secondary fixative after neutral buffered formalin (NBF), providing that exposure to NBF has not exceeded several hours.
5. If used in enclosed tissue processors, it is advisable to follow the manufacturer's instructions regarding flushing and housekeeping procedures.

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