# Shandon Cell-Fixx™ Spray Fixative

**Thermo**  
ELECTRON CORPORATION

Anatomical Pathology USA  
Clinical Diagnostics  
171 Industry Drive  
Pittsburgh, PA 15275, USA  
Tel: 1-800-547-7429  
+1 412 788 1133  
Fax: +1 412 788 1138  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

Anatomical Pathology International  
Clinical Diagnostics  
93-96 Chadwick Road  
Astmoor, Runcorn  
Cheshire, WA7 1PR, UK  
Tel: +44 (0) 1928 566611  
Fax: +44 (0) 1928 565845  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

**Shandon Cell-Fixx™ Spray Fixative** is a convenient and reliable cytology fixative containing Carbowax and alcohol. CELL-FIXX penetrates cells rapidly, preventing air drying. The fixed slides can then be transported, stained or stored.

CELL-FIXX begins fixing immediately. It preserves the cell morphology and nuclear detail. Because it contains alcohol, CELL-FIXX dehydrates the cell, causing it to shrink slightly. This process sharpens nuclear detail and provides the chromatin patterns routinely seen in cytologic preparations.

The Carbowax infiltrates and supports cellular structure, helps maintain nuclear detail and guards against accidental air drying. In addition, when slides prepared with Carbowax are air dried, the cells firmly adhere to the slide preventing loss of cells during the staining process. Slides fixed with CELL-FIXX can be stored for weeks or months with no adverse effect.

#### General Comments:

CELL-FIXX is a cytology fixative that provides rapid fixation of any (GYN or non-GYN) cytological preparation on a microscope slide. Fixed slides are ready for transporting, staining or storage in about one to five minutes depending on room temperature and humidity.

CELL-FIXX, when sprayed on a smear, will fix it rapidly at room temperature, and leave the smear under a waxy protective film that is water and alcohol soluble. This film maintains the integrity and characteristics of the fixed smear until the film is dissolved in the normal staining process.

### **WARNINGS AND PRECAUTIONS FOR IN VITRO DIAGNOSTIC USE**

#### **WARNING - FLAMMABLE LIQUID**

Do not handle or store near heat, sparks, flames or strong oxidants.

Airborne Exposure Limit: 400 ppm

Inhalation of high vapor concentrations may cause headaches, nausea, confusion and drowsiness. Vapors also cause irritation of eyes, nose and throat. Ingestion can cause irritation and necrosis to kidneys and liver. Repeated and prolonged contact of liquid with skin can cause coloring, dryness, and irritation

DISPOSAL should be completed in accordance with federal, state and local regulations

## WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

### Shelf Life

Shandon Cell-Fixx™ has an unopened shelf life of two years from date of manufacture, indicated by the "Use By" date on the label. Once opened, it should remain usable for 6 months or until the "Use by" date if properly stored.

### Storage

Store in areas approved for flammables.

Do not handle or store near heat, sparks, flames, or strong oxidants.

Keep containers closed when not in use.

## INSTRUCTIONS FOR USE

1. Prepare a thin, even smear on a clean microscope slide.
2. **IMMEDIATELY** after preparation, spray a complete and even coating of CELL-FIXX onto the cells from about 10-12 inches away.  
**WARNING – DO NOT LET THE SPECIMEN DRY BEFORE SPRAYING!**
3. Fixed slides will dry in about one to five minutes, at which time they are ready for transport, staining or storage.
4. The wax coating must be removed to ensure proper staining. This is best completed by soaking slides in 95% alcohol for 15 minutes immediately before staining.
5. The sprayed slides may be stored for months with complete preservation.

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