# Shandon Cryomatrix™ Frozen Embedding Medium

**Thermo**  
ELECTRON CORPORATION

Anatomical Pathology USA  
Clinical Diagnostics  
171 Industry Drive  
Pittsburgh, PA 15275, USA  
Tel: 1-800-547-7429  
+1 412 788 1133  
Fax: +1 412 788 1138  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

Anatomical Pathology International  
Clinical Diagnostics  
93-96 Chadwick Road  
Astmoor, Runcorn  
Cheshire, WA7 1PR, UK  
Tel: +44 (0) 1928 566611  
Fax: +44 (0) 1928 565845  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

**Shandon Cryomatrix™** is an outstanding embedding medium used to support tissue for frozen sectioning. It forms a firm bond with the specimen holder that secures tissue during microtomy. Cryomatrix freezes quickly, sections smoothly, and dissolves in water without residue.

## **WARNINGS AND PRECAUTIONS FOR IN VITRO DIAGNOSTIC USE.**

---

---

### **WARNINGS AND PRECAUTIONS**

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions.

**Definitions are at the rear of this booklet**

---

---

### **SHELF LIFE**

Unopened, Shandon Cryomatrix has a shelf life of 4 years. Once opened, it should remain usable for 6 months if kept capped and properly stored.

### **STORAGE**

Shandon Cryomatrix should be stored at room temperature. Keep containers tightly closed when not in use.

### **INSTRUCTIONS FOR USE:**

1. Cover the surface of the specimen holder with a thick layer of Cryomatrix™.
2. Position and orient tissue on holder, and press lightly to eliminate air bubbles beneath the tissue.
3. Apply additional Cryomatrix until the tissue is completely surrounded.
4. Freeze tissue to holder in usual manner until tissue and Cryomatrix appear white and firm.
5. Section specimen in usual manner.
6. No special procedure is required to remove Cryomatrix from sections; it will dissolve during staining.
7. To remove specimen from holder for further processing, dissolved Cryomatrix in water and immediately place specimen in fixative.

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