Shandon EZ-MOUNT™ Mounting Medium
Shandon-Mount™ Mounting Medium
Shandon Synthetic Mounting Medium
Shandon Xylene Substitute Mounting Medium

Thermo
ELECTRON CORPORATION

Anatomical Pathology USA
Clinical Diagnostics
171 Industry Drive
Pittsburgh, PA 15275, USA
Tel: 1-800-547-7429
+1 412 788 1133
Fax: +1 412 788 1138
www.thermo.com/shandon

Anatomical Pathology International
Clinical Diagnostics
93-96 Chadwick Road
Astmoor, Runcorn
Cheshire, WA7 1PR, UK
Tel: +44 (0) 1928 566611
Fax: +44 (0) 1928 565845
www.thermo.com/shandon

SHANDON MOUNTING MEDIUMS

EZ-MOUNT™ is a xylene-based mounting medium for use with glass coverslips. It is clear and colorless and will not discolor with age. The EZ-MOUNT formula contains an anti-oxidant to inhibit stain fading. Its refractive index is 1.495 + 0.005 and it is soluble in both xylene and toluene.

SHANDON-MOUNT™ is a toluene-based mounting medium that can be used as a conventional mounting medium or it can be employed as a liquid coverslip. SHANDON-MOUNT is clear and colorless and will not discolor with age. It contains an anti-oxidant to inhibit stain fading. When used as a liquid coverslip, SHANDON-MOUNT dries overnight for next day filing. SHANDON-MOUNT is soluble in both xylene and toluene.

SHANDON SYNTHETIC MOUNTANT is a toluene-based mounting medium that can be used as a conventional mounting medium. Shandon Synthetic Mountant is clear and colorless and will not discolor with age. It contains an anti-oxidant to inhibit stain fading. Shandon Synthetic Mountant's refractive index is 1.495 + 0.005 and it is soluble in both xylene and toluene.

XYLENE SUBSTITUTE MOUNTANT was specially formulated for its compatibility with Shandon Xylene Substitute, but it will work very well in all routine mounting of specimens. XYLENE SUBSTITUTE MOUNTANT is soluble in both xylene and toluene. Although compatible with Shandon Xylene Substitute, this mounting medium is not soluble in it. See the section on "Removing A Coverslip."

WARNINGS AND PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE

Reasonable care should be taken when using all laboratory reagents. Use with adequate ventilation.

STORAGE:

Store containers in areas approved for flammables. Do not handle or store near heat, sparks, flames or strong oxidants. Keep containers closed when not in use.

WARNING: DO NOT USE XYLENE SUBSTITUTE MOUNTANT WITH PICRIC ACID-BASED STAINS. However, Picric Acid fixation is not an issue.

WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

INSTRUCTIONS FOR USE:

Conventional Coverslipping

Coverslip one slide at a time to prevent brown artifact from forming, and use the minimal amount of medium necessary to effectively bond the coverglass to the slide.

Suggested protocol:

1. Lay the coverslip on an absorbent surface.
2. Place a drop or two of mounting medium toward the edge of the coverslip.
3. Remove a slide from clearing agent.
4. Quickly wipe the back of the slide dry.
5. Invert the slide so that the specimen side faces down.
6. Lower the slide until it reaches the coverslip, starting at one long edge. Gently press until the mounting medium spreads over the entire surface to be joined.
7. Turn over the slide/coverslip assembly.
8. If necessary, drain extra medium by standing the slide on its edge on an absorbent towel.
9. Gently press out any noticeable air bubbles.

NOTE: Thermo has not validated the protocol illustrated here and takes no responsibility for its use. Customers should always validate protocols before placing any reliance on them.

Special Instructions for SHANDON-MOUNT™

SHANDON-MOUNT™ may be employed as a conventional mounting medium or as a liquid coverslip. The following steps describe how to use a liquid coverslip:

1. Drain slide, but leave moist with clearing agent.
2. Place slide on parallel applicator sticks to raise it off a protected level surface.
3. Flow SHANDON-MOUNT onto specimen until it approaches the edges of the slide.
4. Let slide dry naturally on a level surface until tacky (approximately 20 minutes depending on ambient temperature and humidity).
Let stand overnight before filing.

Where improved microscopic resolution is desired, SHANDON-MOUNT may be removed by soaking the slide in xylene. A glass coverslip may then be applied with mounting medium in the usual manner.

Removing a Coverslip

It is sometimes necessary to remove a coverslip. This is usually done so that a slide can be de-stained and re-stained. In this way a slide that was stained improperly can be corrected, or the original stain can be removed to make room for the application of a special stain.

To remove a coverslip do **one** of the following:

1. Soak the slide in xylene until the coverslip falls off.

OR –

2. Place the slide in a 60°C oven for 3 to 4 hours.

OR –

3. Place the slide in a freezer for several minutes.

CAUTION: DO NOT TRY TO REMOVE A COVERSIP BY FORCE.

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