**Shandon Cryochrome™**  
**Frozen Section Embedding Compounds**

**Thermo**  
ELECTRON CORPORATION

Anatomical Pathology USA  
Clinical Diagnostics  
171 Industry Drive  
Pittsburgh, PA 15275, USA  
Tel: 1-800-547-7429  
+1 412 788 1133  
Fax: +1 412 788 1138  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

Anatomical Pathology International  
Clinical Diagnostics  
93-96 Chadwick Road  
Astmoor, Runcorn  
Cheshire, WA7 1PR, UK  
Tel: +44 (0) 1928 566611  
Fax: +44 (0) 1928 565845  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

**SHANDON** frozen section embedding mediums are designed to aid in the identification of tissue samples in the cryostat.

**SHANDON CRYOCHROME™** consists of four bottles of frozen section embedding mediums, colored blue, pink, yellow and green. Tissue on each specimen holder is easily identified by the color of the frozen embedding medium. Color-coded microscope slides in matching colors are also available. Shandon also offers this product in a clear formulation as **SHANDON CRYOMATRIX™**. Also available is **SHANDON M-1** embedding medium, designed to be used for section specimens that cut best at temperatures around -12° to -14° C.

## FOR IN VITRO DIAGNOSTIC USE

### STORAGE

Keep containers closed when not in use.

### SHELF LIFE

Four (4) years. After opening, 6 months if kept tightly capped.

---

## WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

---

**INSTRUCTIONS FOR USE:**

1. Cover the surface of the specimen holder with a thick layer of embedding medium. Use one color for each patient sample, or for each quadrant of a divided specimen.
2. Position and orient tissue on holder and press lightly to eliminate air bubbles beneath the tissue.
3. Apply additional embedding medium until the tissue is surrounded.
4. Freeze tissue to holder in usual manner until tissue and medium are firm. The color of the medium will lighten considerably but will still be easily identifiable.
5. Section in usual manner. If desired, place sections on microscope slides of the color matching the medium used for that specimen.
6. No special procedure is required to remove embedding medium from sections; it will dissolve in alcohol or aqueous solutions during the staining procedure.
7. To remove specimen from holder for further processing, dissolve medium in water and immediately place specimen in fixative.

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