# Shandon Rapid-Chrome™ Frozen Section Staining Kit

**Thermo**  
ELECTRON CORPORATION

Anatomical Pathology USA  
Clinical Diagnostics  
171 Industry Drive  
Pittsburgh, PA 15275, USA  
Tel: 1-800-547-7429  
+1 412 788 1133  
Fax: +1 412 788 1138  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

Anatomical Pathology International  
Clinical Diagnostics  
93-96 Chadwick Road  
Astmoor, Runcorn  
Cheshire, WA7 1PR, UK  
Tel: +44 (0) 1928 566611  
Fax: +44 (0) 1928 565845  
[www.thermo.com/shandon](http://www.thermo.com/shandon)

## RAPID-CHROME™ FROZEN SECTION STAINING KIT

This special staining kit simplifies the routine staining of frozen tissue sections from a cryostat.

With the Rapid-Chrome Staining Kit, all the reagents for permanently stained frozen sections are supplied in one convenient carrying tray. The die-cut slots in the tray base accommodate each jar lid to reduce the problem of cross-contamination of reagents. In addition, the reagents are sequentially aligned so that the technician can dip the microscope slide progressively from jar to jar.

This kit can stain approximately 100 slides, although the actual number may vary due to the type and thickness of the specimen.

### WARNINGS AND PRECAUTIONS FOR IN VITRO DIAGNOSTIC USE

**WARNING: CONTAINS FLAMMABLE REAGENTS**

### STORAGE

Keep away from direct sunlight. Store at room temperature. Do not handle or store near heat, sparks, flames or strong oxidants.

Lot number and expiration date are printed on the tray, not on the individual components

---

### WARNINGS AND PRECAUTIONS

See Material Safety Data Sheets for warnings and precautions, as well as R & S code definitions. The following symbols are used on the product labels and packaging:

Definitions are at the rear of this booklet

---

**NOTE:** This kit contains Reagent Alcohol (100% and 95%) and Xylene, which are considered general laboratory reagents.

## INSTRUCTIONS FOR USE

### A. Kit Preparation

1. All solutions required are provided in the tray and are ready to use full strength.
2. If less intense cytoplasmic staining is desired, Eosin –Y can be diluted with 80 – 85% alcohol (ethanol, reagent, or methanol).
3. Bluing Reagent should be changed as often as necessary to retain its alkalinity.
4. Other solutions requiring periodic replacement are:
  - Distilled Water – frequently
  - Alcohols – as needed (ethanol or reagent)
  - Xylene – as needed

### B. Specimen Preparation

1. To minimize the formation of artifacts and distortion due to air-drying and autolysis, freeze tissue as quickly as possible.
2. Adhere the tissue to the specimen holder ("microtome chuck") using Shandon Cryomatrix™.  
The use of Cryomatrix increases the embedding mass and tends to reduce any distortions occurring during cutting and handling.
3. Cut and mount the frozen section in the usual manner.
4. Place slide in Rapid-Fixx™ immediately.

### C. Staining Procedure

1. Remove all lids and place each in its slot to the right of the jar.
2. Place slides in Rapid-Fixx™ without delay.
3. Process slides through solutions following the directions indicated by the arrows.
4. A suggested staining protocol can be found on the following page.\* Staining times may be adjusted to suit individual preferences in stain intensity.
5. Coverslip the slide with Shandon Mounting Medium.
6. Replace distilled waters, alcohols, xylene and Bluing Reagent if necessary.
7. Replace all lids on their original jars to prevent possible contamination.

\* NOTE: Thermo has not validated the protocol illustrated here and takes no responsibility for its use. Customers should always validate protocols before placing any reliance on them.

FROZEN SECTION STAINING PROTOCOL		
Step	Solution	Time
1	Rapid-Fixx™	5 – 7 seconds
2	Distilled Water	5 – 10 dips
3	Hematoxylin	1 minute
4	Distilled Water	5 – 10 dips
5	Bluing Reagent	3 dips
6	95% Alcohol	5 – 7 dips
7	Eosin-Y	15 seconds
8	95% Alcohol	5 – 7 dips
9	100% Alcohol	5 – 7 dips
10	100% Alcohol	5 – 7 dips
11	Xylene	5 – 7 dips
12	Xylene	5 – 7 dips

**ACTIVE INGREDIENTS**

## 100% Alcohol:

Specially Denatured Alcohol Formula 3.A 95 Parts by volume  
Isopropyl Alcohol 5 Parts by volume

**DANGER. POISON. VAPOR HARMFUL. MAY BE FAATAL OR CAUSE BLINDNESS IF SWALLOWED. CANNOT BE MADE NON-POISONOUS. FLAMMABLE.**

## 95% Alcohol:

Reagent alcohol

**DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED. CANNOT BE MADE NON-POISONOUS. FLAMMABLE.**

## Rapid-Fixx™

Methanol 75.0%  
Formaldehyde 20.0%  
Glacial Acetic Acid 5.0%

**WARNING: STRONG SENSITIZER. AVOID PROLONGED OR REPEATED EXPOSURE. DANGER. POISON. VAPOR HARMFUL. MAY BE FATAL OR CAUSE BLINDNESS IF SWALLOWED. CANNOT BE MADE NON-POISONOUS. FLAMMABLE.**

## Gill 3 Hematoxylin:

Aluminum sulfate 3.0%  
Hematoxylin 0.68%

**WARNING: HARMFUL OR FATAL IF SWALLOWED.**

## Bluing Reagent:

Sodium bicarbonate 0.2%  
Lithium carbonate 0.04%

**DANGER. POISON. VAPOR HARMFUL. MAY BE FAATAL OR CAUSE BLINDNESS IF SWALLOWED. CANNOT BE MADE NON-POISONOUS. FLAMMABLE.**

Eosin-Y:

Eosin-Y 0.5%

**DANGER. POISON. VAPOR HARMFUL. MAY BE FAATAL OR CAUSE BLINDNESS IF SWALLOWED. CANNOT BE MADE NON-POISONOUS. FLAMMABLE.**

## Xylene:

Xylene 100%

**DANGER. HARMFUL OR FATAL IF SWALLOWED. IF SWALLOWED, DO NOT INDUCE VOMITING. CALL A PHYSICIAN IMMEDIATELY. VAPOR HARMFUL. FLAMMABLE.**

## Shandon Mounting Medium:

Methacrylate polymer 57%

**DANGER. HARMFUL OR FATAL IF SWALLOWED. IF SWALLOWED, DO NOT INDUCE VOMITING. CALL A PHYSICIAN IMMEDIATELY. VAPOR HARMFUL. FLAMMABLE.**

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