

Electron Capture Dissociation on the Finnigan LTQ FT–Preserving Post-translational Modifications during Peptide Fragmentation

Robert Malek, Wolfgang Metelmann-Strupat, Martin Zeller and Helmut Muenster, Thermo Electron Corporation, Bremen, Germany

Key Words

- Finnigan™ LTQ FT™
- ECD
- FTICR-MS
- Post-translational Modification (PTM)
- Proteomics

Introduction

Electron capture dissociation (ECD) has recently evolved as an alternate activation method, especially for peptide and protein sequencing with Fourier-transform ion cyclotron resonance-mass spectrometry (FTICR-MS). With ECD, multiply charged cations are irradiated with low energy electrons produced by an emitter cathode behind the ICR cell. Electron capture produces a radical cation $[M+nH]^{(n-1)+\bullet}$ which can dissociate by a rapid, facile fragmentation of the N-C α bond of the peptide chain, producing mainly c- and z*-type fragment ions^[1]

In contrast to ECD, the most commonly used activation methods such as collision induced dissociation (CID) and infrared multiphoton dissociation (IRMPD) induce dissociation by vibrational excitation of the precursor. Since these activation methods induce “ergodic” processes, i.e. they add internal energy to the precursor slower than the rate of energy randomization, usually cleavage of the weakest bonds within the precursor is observed. Within peptides, the backbone amide bond has the lowest energy barrier to dissociation and predominantly b- and y-type fragment ions are formed. However, substituents added in co- and post-translational modifications often have lower energy barriers than those of backbone cleavage. This can result in more complex tandem mass spectra and potentially the loss of the information on the attachment site of these substituents.

Due to the non-ergodic nature of ECD, co- and post-translational modifications are preserved. ECD allows site specific analysis of phosphorylation,^[2,3] O- and N-linked glycosylation,^[4-7] and sulfation^[8] ECD holds much promise as a supplementary dissociation technique to CID for unambiguous protein identification, *de novo* sequencing^[9] and detailed protein characterization.^[10]


ECD shows higher degrees of fragmentation, allowing the distinction of leucine and isoleucine residues^[11-13] and even between D- and L-amino acids!^[14] ECD can also provide additional information of the amino acid composition by careful examination of side chain losses!^[15-17]

However, the overall efficiency of ECD is typically lower than that obtained with CAD. Long ion accumulation, activation, and detection times in the ICR cell together with the need for the addition of multiple spectra makes the use of ECD for on-line separation of complex peptide or protein mixtures quite challenging.

The Finnigan LTQ FT hybrid linear ion trap/FTICR mass spectrometer provides a means to perform ECD analysis on complex mixtures of peptides and proteins on an LC timescale. The use of ion accumulation and isolation in the linear ion trap prior to transfer to the ICR cell for ECD and product ion detection allows improved cycle times (< 1.5 seconds), resulting in the acquisition of multiple ECD scans for each peak eluted. High cleavage efficiency is achieved by the use of a cathode with large electron emitting surface^[18].

The cathode is positioned at the back of the ICR cell, slightly off-axis to allow the simultaneous use of ECD and IRMPD (see below).

Figure 1 shows a single scan ECD MS/MS spectrum of the doubly charged precursor of substance P at m/z 674.37. The spectrum exhibits intense ECD fragment



ion peaks. Bearing in mind that the cyclic structure of proline does not allow formation of c- and z-type fragments, all possible N-C α bonds are cleaved, allowing even *de novo* sequencing of peptides with unknown amino acid sequences.

The CID MS/MS spectrum of substance P is shown in Figure 2. The doubly charged peptide precursor ions were subjected to CID in the linear ion trap and the fragment ions were transferred into the ICR cell and detected. The spectrum looks somewhat more complex compared to the ECD spectrum. The fragment ion peaks of this spectrum

are sufficient to identify substance P in a database search, but *de novo* sequencing would be a challenge.

As mentioned above, labile post-translational modifications are preserved in ECD. This is due to the fact that the amide carbonyl has a higher H $^+$ affinity than the amino acid side chains and, therefore, the backbone N-C α bonds are cleaved preferentially.

Figure 3 shows the ECD MS/MS spectrum of an O-glycosylated peptide carrying one N-acetylglucosamine (GlcNAc) moiety attached to one of its serines.




Figure 1: ECD spectrum of substance P, single scan at 100,000 resolution (FWHM)




Figure 2: CID spectrum of substance P, single scan at 100,000 resolution (FWHM)

Interestingly, ECD of this glycosylated peptide gives rise to z- and z[•]-type fragment ion peaks in addition to one c-type fragment ion peak. The fragment ion peaks z₄ to z₇ are one Dalton heavier than expected and can be explained by direct H[•] capture or H rearrangement and is a previously observed minor pathway in ECD!¹⁹

However, the mass difference between the z₆ and z₇ fragment ion peaks corresponds to the mass of a serine carrying a GlcNAc moiety, whereas the mass difference

between the z₄ and z₅-fragment ion peaks corresponds to an unmodified serine residue. Although no peak for z₁₁[•]- fragment ions can be seen in the spectrum, the mass difference between the z₁₀[•]- fragment ion peak and the mass of the intact peptide corresponds to unmodified serine and proline.

The CID spectrum of the doubly charged precursor of this glycosylated peptide is shown in Figure 4.




Figure 3: ECD spectrum of the O-glycosylated peptide




Figure 4: CID spectrum of the O-glycosylated peptide

It can be clearly seen in the spectrum that the major fragmentation pathway under CID conditions is the loss of the neutral sugar moiety. The spectrum also reveals fragment ions with the sugar moiety still attached and, in addition, fragments with the loss of the neutral sugar. The abundant peaks for y_0 -fragment ions can be explained by the presence of proline in the peptide sequence which is known to be prone to facile fragmentation under CID conditions.

Conclusion

ECD is an effective alternate dissociation method, especially for the analysis of post translational modifications on an LC time scale. With the LTQ FT, additional confidence in the results is provided by accurate mass assignment of the fragment ions.

References:

- [1] Zubarev, R.A., Kelleher, N.L., and McLafferty, F.W., Electron capture dissociation of multiply charged protein cations. A nonergodic process, *J. Am. Chem. Soc.* 1998, 120, 3265-3266.
- [2] Stensballe, A., Jensen, O.N., Olsen, J.V., Haselmann, K.F., and Zubarev, R.A., Electron capture dissociation of singly and multiply phosphorylated peptides, *Rapid Commun. Mass Spectrom.* 2000, 14, 1793-1800.
- [3] Shi, S.D., Hemling, M.E., Carr, S.A., Horn, D.M., Lindh, I., and McLafferty, F.W., Phosphopeptide/phosphoprotein mapping by electron capture dissociation mass spectrometry, *Anal. Chem.* 2001, 73, 19-22.
- [4] Mirgorodskaya, E., Roepstorff, P., and Zubarev, R.A., Localization of O-glycosylation sites in peptides by electron capture dissociation in a Fourier transform mass spectrometer, *Anal. Chem.* 1999, 71, 4431-4436.
- [5] Hakansson, K., Cooper, H.J., Emmett, M.R., Costello, C.E., Marshall, A.G., and Nilsson, C.L., Electron capture dissociation and infrared multiphoton dissociation MS/MS of an N-glycosylated tryptic peptide to yield complementary sequence information, *Anal. Chem.* 2001, 73, 4530-4536.
- [6] Budnik, B.A., Haselmann, K.F., Elkin, Y.N., Gorbach, V.I., and Zubarev, R.A., Applications of electron-ion dissociation reactions for analysis of polycationic chitoooligosaccharides in fourier transform mass spectrometry, *Anal. Chem.* 2003, 75, 5994-6001.
- [7] Mormann, M., Macek, B., de Peredo, A.G., Hofsteenge, J., and Peter-Katalinic, J., Structural studies on protein O-fucosylation by electron capture dissociation, *Int. J. Mass Spectrom.* 2004, 234, 11-21.
- [8] Kelleher, N.L., Zubarev, R.A., Bush, K., Furie, B., Furie, B.C., McLafferty, F.W., and Walsh, C.T., Localization of labile posttranslational modifications by electron capture dissociation: the case of gamma-carboxyglutamic acid, *Anal. Chem.* 1999, 71, 4250-4253.
- [9] Budnik, B.A., Olsen, J.V., Egorov, T.A., Anisimova, V.E., Galkina, T.G., Musolyamov, A.K., Grishin, E.V., and Zubarev, R.A., De novo sequencing of antimicrobial peptides isolated from the venom glands of the wolf spider *Lycosa singoriensis*, *J. Mass Spectrom.* 2004, 39, 193-201.
- [10] Kjeldsen, F., Haselmann, K.F., Budnik, B.A., Sorensen, E.S., and Zubarev, R.A., Complete characterization of posttranslational modification sites in the bovine milk protein PP3 by tandem mass spectrometry with electron capture dissociation as the last stage, *Anal. Chem.* 2003, 75, 2355-2361.
- [11] Kjeldsen, F., Haselmann, K.F., Budnik, B.A., Jensen, F., and Zubarev, R. A., Dissociative capture of hot (3-13 eV) electrons by polypeptide polycations: an efficient process accompanied by secondary fragmentation, *Chem. Phys. Lett.* 2002, 356, 201-206.
- [12] Kjeldsen, F., Haselmann, K.F., Sorensen, E.S., and Zubarev, R.A., Distinguishing of Ile/Leu amino acid residues in the PP3 protein by (hot) electron capture dissociation in Fourier transform ion cyclotron resonance mass spectrometry, *Anal. Chem.* 2003, 75, 1267-1274.
- [13] Kjeldsen, F., and Zubarev, R., Secondary losses via gamma-lactam formation in hot electron capture dissociation: A missing link to complete de novo sequencing of proteins?, *J. Am. Chem. Soc.* 2003, 125, 6628-6629.
- [14] Adams, C.M., Kjeldsen, F., Zubarev, R.A., Budnik, B.A., and Haselmann, K.F., Electron capture dissociation distinguishes a single D-amino acid in a protein and probes the tertiary structure, *J. Am. Soc. Mass Spectrom.* 2004, 15, 1087-1098.
- [15] Cooper, H.J., Hudgins, R.R., Hakansson, K., and Marshall, A.G., Characterization of amino acid side chain losses in electron capture dissociation., *J. Am. Soc. Mass Spectrom.* 2002, 13, 241-249.
- [16] Haselmann, K.F., Budnik, B.A., Kjeldsen, F., Polfer, N.C., and Zubarev, R.A., Can the (Mo - X) region in electron capture dissociation provide reliable information on amino acid composition of polypeptides?, *Eur. J. Mass Spectrom.* 2002, 8, 461-469.
- [17] Cooper, H.J., Hakansson, K., Marshall, A.G., Hudgins, R.R., Haselmann, K.F., Kjeldsen, F., Budnik, B.A., Polfer, N.C., and Zubarev, R.A., Letter: The diagnostic value of amino acid side-chain losses in electron capture dissociation of polypeptides. Comment on: "Can the (M-X) region in electron capture dissociation provide reliable information on amino acid composition of polypeptides?," *Eur. J. Mass Spectrom.* 8, 461-469 (2002), *Eur. J. Mass Spectrom.* 2003, 9, 221-222.
- [18] Tsybin, Y.O., Hakansson, P., Budnik, B.A., Haselmann, K.F., Kjeldsen, F., Gorshkov, M. and Zubarev, R.A., Improved low-energy electron injection systems for high rate electron capture dissociation in Fourier transform ion cyclotron resonance mass spectrometry, *Rapid Commun. Mass Spectrom.* 2001, 15, 1849-1854.
- [19] Zubarev, R.A., Reactions of polypeptides ions with electrons in the gas phase, *Mass Spectrom. Rev.* 2003, 22, 57-77.

In addition to these offices, Thermo Electron Corporation maintains a network of representative organizations throughout the world.

Australia
+61 2 8844 9500

Austria
+43 1 333 50340

Belgium
+32 2 482 30 30

Canada
+1 800 532 4752

China
+86 10 5850 3588

France
+33 1 60 92 48 00

Germany
+49 6103 4080

India
+91 22 2778 1101

Italy
+39 02 950 591

Japan
+81 45 453 9100

Latin America
+1 512 251 1503

Netherlands
+31 76 587 98 88

Scandinavia
+46 8 556 468 00

South Africa
+27 11 570 1840

Spain
+34 91 657 4930

Switzerland
+41 61 48784 00

UK
+44 1442 233555

USA
+1 800 532 4752

www.thermo.com

©2005 Thermo Electron Corporation. All rights reserved. ProSight PTM is a trademark of the University of Illinois at Urbana-Champaign. All other trademarks are the property of Thermo Electron Corporation and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Electron (Bremen) GmbH is certified DIN EN ISO 9001:2000

AN61885_E 10/05S

Thermo
ELECTRON CORPORATION