

# Isolation of RNA Using Thermo Scientific Rotors and Sorvall Discovery SE Microultracentrifuges

Stephanie Noles, Ph.D., Thermo Fisher Scientific, Asheville, NC

## KEY WORDS

- RNA Isolation
- Cesium Chloride
- Guanidine Thiocyanate
- S150-AT Fixed-Angle Rotor
- S55-S Swinging Bucket Rotor
- Microultracentrifuge
- Microcentrifuge

## Introduction

Many protocols exist for the extraction and purification of total RNA from cells or tissue. The key to obtaining high quality RNA is the rapid lysis of the cells, the thorough denaturation of endogenous nucleases, and ultimately, the efficient removal of these nucleases. Early methods of RNA purification relied on the use of guanidinium thiocyanate (GTC) for the rapid denaturation of cellular proteins and extraction of RNA by ultracentrifugation through a cesium chloride (CsCl) cushion or a single step procedure that uses a mixture of GTC and phenol-chloroform. Recent techniques and commercially available kits have been developed around these methods that allow for a more rapid process for RNA isolation. Using microcentrifuges, these methods have been found to decrease the time for RNA isolation from a few hours to less than an hour. This brief will provide protocols for the traditional methods of RNA isolation using fixed-angle and swinging bucket rotors in a Thermo Scientific Sorvall Discovery SE microultracentrifuge, and will also identify commercially available kits that are compatible with Thermo Scientific microcentrifuges.

## Procedures

### PROTOCOL 1: Rapid Isolation of RNA from Cultured Cells Using the Thermo Scientific S55-S Swinging Bucket Rotor and Sorvall® Discovery™ SE Microultracentrifuge<sup>5</sup>

This report describes the use of GTC and centrifugation through a CsCl cushion with a S55-S swinging bucket rotor to isolate RNA from

cells grown in culture. This procedure has been developed for working with small samples and offers a high yield (>90%) of RNA in a short period of time (approximately 3 hours).

## Lysis

1. Remove media from cell culture plate.
2. Add 600-1000  $\mu\text{L}$  of GTC solution (4.0 M guanidine thiocyanate, 0.03 M NaOAc, 1%  $\beta$ -mercaptoethanol) to each well of the culture plate to lyse.
3. Lift cells with a cell scraper and homogenize the slurry by aspirating 4-5 times through a 22-gauge needle.

## Isolation

1. Pre-load 2.2 mL polyallomer tubes (Cat. No. 45240) with 550  $\mu\text{L}$  of 5.7 M CsCl and mark solution level with a waterproof pen (Figure 1).
2. Layer samples onto the CsCl cushion and fill tubes to within 3 mm of the top of the tube with GTC solution.
3. Place tubes in an S55-S swinging bucket rotor and centrifuge in a Sorvall Discovery SE microultracentrifuge at 55,000 rpm (258,826  $\times$  g) for 3.0 – 3.2 hours at 22 °C. (Note: Non-precipitating solutions up to 1.7 g/mL can be run in the S55-S rotor without a reduction in rotor speed. Consult rotor manual for instructions on rotor de-ration and how to reduce the chance of CsCl precipitation).
4. Following centrifugation, use a sterile Pasteur pipette to remove the top layer of solution down to the interface reference mark (Figure 2).



Thermo Scientific Sorvall Discovery SE Microultracentrifuge




Thermo Scientific S55-S Swinging Bucket Rotor

5. Using a second sterile pipette, remove remaining solution down to approximately 250  $\mu\text{L}$ .
6. Carefully, cut off the top of each tube just below the reference mark, pour off the remaining solution, and blot the tubes dry on a clean paper towel.

## Recovery

1. Rinse RNA pellet  $\times$  2 with 75  $\mu\text{L}$  of 70% isopropanol/0.3 M sodium acetate and blot the tubes dry on a clean paper towel.
2. Add 100  $\mu\text{L}$  of DEPC-treated water to each tube and freeze (-70 °C) for a minimum of 45-60 minutes. (Note: DEPC is a suspected carcinogen; follow manufacturers recommended handling procedures).

- Resuspend RNA by pipetting up and down several times.
- Transfer solution to a 1.5 mL microcentrifuge tube.
- Wash the ultracentrifuge tube with a second 100  $\mu$ L aliquot of 0.25 volumes of 7.5 M ammonium acetate and 2.6 volumes of 100% ethanol.
- Collect precipitated RNA by centrifugation at  $>10,000 \times g$  for 20 minutes at 4 °C in a Thermo Scientific Heraeus Pico or Fresco refrigerated microcentrifuge.
- Wash the pellet with cold 70% ethanol and vacuum dry.
- Resuspend RNA pellet in a small volume of DEPC-treated water.



**Figure 1.** Polyallomer tube containing the sample (GIT solution) on top of the CsCl cushion.


due to crystallization during centrifugation when 5.7M CsCl is used.)

### Conclusion

The isolation and purification of RNA can be accomplished using multiple methods. This brief describes traditional approaches for RNA purification using ultracentrifugation with swinging bucket and fixed-angle rotors. A swinging bucket rotor is often used because the RNA is collected at the bottom of the tube. Although centrifugation in a vertical or fixed-angle rotor requires less time for RNA purification, it may result in materials collecting or coming in contact with the sides of the tube, thus potentially contaminating the RNA. Alternatively, commercially available kits (ie, Promega® Ambion®, Qiagen™) have been developed that are compatible with Thermo Scientific microcentrifuges and may offer enhanced time savings, while increasing RNA yield and purity.

### References

- Chirgwin, J.J., A.E. Przbyla, R.J. MacDonald, and W.J. Rutter. 1979. Isolation of



**Figure 2.** Withdrawal of the supernatant following centrifugation. The figure illustrates the depth of the pipette tips and the cut position relative to the reference mark placed on the tube delineating the interface between the sample and the CsCl cushion.

Biologically Active Ribonucleic Acid From Sources Enriched in Ribonuclease. *Biochemistry* 18: 5294.

- Glisin, V., R. Crkvenjakov, and C. Byrus. 1973. Ribonucleic Acid Isolated by Cesium Chloride Centrifugation. *Biochemistry* 13:2633.
- Norton, J. Rapid Isolation of RNA from Cultured Cells Using the S55-S Swinging Bucket Rotor and Discovery Microultracentrifuge. Application Brief, S00266.
- Sambrook J., et al.: *Molecular Cloning* (2nd edition); Cold Spring Harbor Laboratory (1989).
- Separation of total RNAs from a homogenated cell line with the S150-AT fixed angle rotor in less than 2 hours. Diverse Applications for Discovery M120 & M150 Micro-ultracentrifuges. Application Brief, S00242.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

#### North America:

**USA / Canada**  
+1 866 984 3766

#### Europe:

**Austria**  
+43 1 801 40 0

**Belgium**  
+32 2 482 30 30

**France**  
+33 2 2803 2000

**Germany national toll free**  
08001-536 376

**Germany international**  
+49 6184 90 6940

**Italy**  
+39 02 02 95059 341

**Netherlands**  
+31 76 571 4440

**Nordic countries**  
+358 9 329 100

**Russia / CIS**  
+7 (812) 703 42 15

**Spain/Portugal**  
+34 93 223 09 18

**Switzerland**  
+41 44 454 12 12

**UK/Ireland**  
+44 870 609 9203

#### Asia:

**China**  
+86 21 6865 4588 or +86 10 8419 3588

**India**  
+91 22 6716 2200

**Japan**  
+81 45 453 9220

**Other Asian countries**  
+852 2885 4613

**Countries not listed**  
+49 6184 90 6940 or +33 2 2803 2000

[www.thermo.com/centrifuge](http://www.thermo.com/centrifuge)

© 2008 Thermo Fisher Scientific Inc. All rights reserved. Re-Seal is a trademark of Seton Scientific. Promega is a registered trademark of Promega Corporation. Ambion is a registered trademark of Ambion, Inc. Qiagen is a trademark of Qiagen. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.