

High-Throughput Qualitative and Quantitative Analysis of High Resolution Accurate Mass Data for Early Drug Discovery

Lester Taylor¹, Karen Salomon¹, Caroline Ding¹, Hans Pfaff², David Wright¹, Jim Shofstahl¹, Prakash Chand¹ and Zesmaa Morimoto¹

1. Thermo Fisher Scientific, San Jose, CA, USA

2. Thermo Fisher Scientific, Bremen, Germany



Overview

Purpose: We describe the use of instrument software for metabolism screening which controls acquisition and processing of ultra-high resolution accurate mass data from the Thermo Scientific Exactive LC-MS instrument. **Methods:** The software incorporates several novel algorithms for data processing, including advanced component detection, parameter-less peak integration and intelligent elemental composition determination. **Results:** The software can be used to screen and report the nature and extent of parent drug metabolism in an automated fashion. At the same time, information is obtained for putative metabolites and these results are also output.

Introduction

The availability of full scan ultra high resolution LC-MS instruments allow the analysis of samples for drug discovery, metabolism and pharmacokinetics (PK) without the need to set up specific methods for target compounds. After acquisition, significant data processing is required to mine the data for parent drug and putative metabolites together with an assessment of quantitative and qualitative changes for time course experiments such as used in *in vitro* and *in vivo* compound assessment studies. We describe the use of software which will automatically process and directly report summary PK data.

Methods

In order to automatically process full scan ultra high resolution data it is necessary to employ software algorithms which will specifically manipulate data in order to identify compounds of interest at mass-to-charge values which have not necessarily been pre-defined. This allows the detection and reporting of qualitative and quantitative information for drug metabolites in time course studies in which a parent drug is specified. Components of the overall data processing methodology include: advanced background subtraction, and peak picking using parameter-less peak detection. Sophisticated elemental composition determination is used to identify parent drug related metabolites.

Workflow


Prototype software has been developed to acquire and process high resolution, accurate mass metabolite data. A schematic of the workflow is shown in Figure 1. After the raw data files are acquired, a novel algorithm is used to detect all components having an acceptable chromatographic profile and a master peak list is created of all components detected. Peak area calculations are performed on these components using a parameter-less peak detection algorithm. Another algorithm is used to determine which components are isotopes or adducts of other components and these peaks are labeled accordingly in the master peak list.

The next step is to compare the peak list of one set of time points to a reference peak list – either a time zero data file or a control file. Any peaks common to both lists are labeled as such and are most likely not metabolites.

Metabolite finding is carried out on mass peaks in a sample that are different from the reference file. Metabolite finding will look at each un-annotated entry in the master peak list for known biotransformations based on arithmetic difference. Peaks matching an expected modification will be labeled with that modification in the master peak list. Peaks not matching an expected modification remain un-annotated, but are not removed from the master peak list.

Plots of drug/metabolite concentration as a function of time are generated using the annotated peak lists for each time course sample. Relative quantitation of target compounds (parent drug and/or metabolites) is based on a user-defined reference. Confirmation of putative metabolites is provided with an elemental composition calculation.

FIGURE 1. High-throughput quan/qual workflow



Automated Component Detection

The algorithm for advanced component detection is designed to remove noise and baseline signals from a data file. It can detect and separate components from each other even if their concentrations are different and their retention times are similar. Processing can be applied to chromatographic full scan high resolution MS spectra, i.e. LC-MS data files. There is a requirement that a minimum of three (or more) MS spectra should be acquired for each eluting component.

The algorithm performs a scan-by-scan analysis of the data file and searches for mass peaks within each scan. For each mass peak present, a mass chromatogram is created. The mass chromatogram of each peak within the scan is then evaluated to determine whether it has an acceptable chromatographic profile (Figure 2). For those masses with an acceptable profile, the centroid retention time is calculated and the mass is added to an event list with the centroid retention time. (The centroid retention time is used because it is more accurate.) The analysis is carried out for each scan in the data file to create a comprehensive event list of all components.

FIGURE 2. Examination of mass peaks in one scan at RT = 90.17 minutes. Mass peaks with acceptable chromatographic profiles are labeled with blue circles and mass peaks with unacceptable chromatographic profiles are labeled with red crosses

