

Certificate of Analysis

Product Information:

Orion Number	080113
Description	PROBE SERVICE KIT FOR 810/850 PROBES
Lot Number	HQ1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	1.28 +/- 0.02 M KCL

Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	1.28	Moles/Liter Potassium Chloride

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 10/28/03

Approved By:



Quality Control Supervisor