

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	011008
<b>Description</b>	COND STD 100 US 5X60ML
<b>Lot Number</b>	HQ1
<b>Expiration Date</b>	10/05
<b>Method of Analysis</b>	Conductivity
<b>Acceptance Range</b>	100.0 uS/cm +/- 2.0 uS/cm @ 25.0

## Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	100.2	uS/cm at 25 deg. C

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Chloride	2202	NA

**Certificate Date:** 11/10/03

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**