

Certificate of Analysis

Product Information:

Orion Number	940909
Description	TISAB II (ONE GALLON)
Lot Number	HO1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	pH 5.0 to 5.5/Fluoride < 0.10 ppm

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.030	ppm Fluoride
2	5.12	pH

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Dihydrogen Phosphate	186 I	f
Disodium Hydrogen Phosphate	186 II	f
Potassium Hydrogen Phthalate	185	g

Certificate Date: 12/12/03

Approved By:

Edward Henning

Quality Control Supervisor