

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	011007
<b>Description</b>	COND/TDS STD 1413us 5X60 ML
<b>Lot Number</b>	HS1
<b>Expiration Date</b>	08/05
<b>Method of Analysis</b>	Conductivity
<b>Acceptance Range</b>	1413 uS/cm +/- 7 uS/cm @ 25.0 C

## Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	1275	uS/cm at 20 deg. C

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Chloride	2202	NA

**Certificate Date:** 12/17/03

**Approved By:**



**Quality Control Supervisor**