

Certificate of Analysis

Product Information:

| | |
|---------------------------|--------------------------|
| Orion Number | 911025 |
| Description | PH 10 BUFFER IN 25 PACK. |
| Lot Number | HR1C |
| Expiration Date | 09/05 |
| Method of Analysis | Potentiometric |
| Acceptance Range | 10.01 +/- 0.02 pH UNITS |

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------|
| 1 | 10.018 | pH |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|--------------------------------|--------------|---------------|
| Disodium Hydrogen Phosphate | 186 II | f |
| Potassium Dihydrogen Phosphate | 186 I | f |
| Sodium Bicarbonate | 191 | b |
| Sodium Carbonate | 192 | b |

Certificate Date: 12/23/03

Approved By:



Quality Control Supervisor