

Certificate of Analysis

Product Information:

Orion Number	011006
Description	COND/TDS STD 12.9MS/CM 5x60ML
Lot Number	H01
Expiration Date	12/05
Method of Analysis	Conductivity
Acceptance Range	0.1000M+/-0.0005M

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.09977	Moles/Liter Potassium Chloride
2	12.87	mS/cm at 25 deg. C

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Chloride	2202	NA

Certificate Date: 12/23/03

Approved By:

Edward Henning

Quality Control Supervisor