

Certificate of Analysis

Product Information:

Orion Number	950202
Description	CO2 INTERNAL FILL SOLN
Lot Number	HP1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	0.10 M +/- 0.01 M KCL

Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.10	Moles/Liter Potassium Chloride

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 11/11/03

Approved By:



Quality Control Supervisor