

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	080013
<b>Description</b>	PROBE SERVICE KIT Electrolyte Solution
<b>Lot Number</b>	HP1
<b>Expiration Date</b>	11/05
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	9.80 to 10.10 pH units @ 25.0 C

## Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	9.88	pH

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	f
Potassium Dihydrogen Phosphate	186 I	f
Potassium Hydrogen Phthalate	185	g

**Certificate Date:** 11/14/03

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**