

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	900062
<b>Description</b>	OPTIMUM RESULTS B RFS 5 X 60ML
<b>Lot Number</b>	FW1
<b>Expiration Date</b>	04/30/03
<b>Method of Analysis</b>	Potentiometric

## Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	0.020	Moles/Liter Potassium Chloride

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

**Certificate Date:** 04/26/01

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**