

Certificate of Analysis

Product Information:

Orion Number	9104CB
Description	PH 4 BUFFER, CUBITAINER 19L
Lot Number	FW1
Expiration Date	04/30/03
Method of Analysis	Potentiometric

Actual Analysis:


<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	4.015	pH

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	f
Potassium Dihydrogen Phosphate	186 I	f
Potassium Hydrogen Phthalate	185	g

Certificate Date: 04/24/01

Approved By:



Quality Control Supervisor