

Certificate of Analysis

Product Information:

Orion Number 080314
Description ELECTROLYTE, D.O. PROBE 50ML BOTTLE
Lot Number 0172-2
Expiration Date 2 Years after Receipt of Product
Method of Analysis Conductivity

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------|
| 1 | 85 +/- 5 | mS/cm at 25 deg. C |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|--------------------------|--------------|---------------|
| Potassium Chloride | 2202 | NA |

Certificate Date: 08/08/01

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245967-001 Rev. B

Thermo
ELECTRON CORPORATION