

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/03/2004

Reviewed on 02/03/2004

## 1 Identification of substance

- **Product details**
- **Trade name:** DPD No.1
- **Product Use:** Reagent
- **Part of Orion:**  
AC2035  
AC2070  
AC2071  
AC2099
- **MSDS p/n:** 233236-001 Revision E
- **Manufacturer/Supplier:**  
Thermo Electron Corporation  
Environmental Instruments  
Water Analysis Group  
166 Cummings Center  
Beverly, MA 01915-6110  
USA  
Phone: 978-232-6000  
Made in Germany

## 2 Composition / Data on components

- **Chemical characterization**
- **Description:** Mixture of the substances listed below with nonhazardous additions.

· CAS No.	Description	Xn; R 20/21/22-36	0.1-1.0%
CAS: 6283-63-2 EINECS: 228-500-6	N,N-diethylbenzene-1,4-diammonium sulphate	<input checked="" type="checkbox"/>	

- **Additional information:**  
For the wording of the listed risk phrases refer to section 16.  
The remaining 99% contains no hazardous ingredients according to OSHA definition.

## 3 Hazards identification

- **Hazard description:** Not applicable.
- **Information pertaining to particular dangers for man and environment:** Not applicable.
- **Classification system:**  
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.
- **NFPA ratings (scale 0 - 4)**  
Health = 0  
Fire = 0  
Reactivity = 0

## 4 First aid measures

- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:** Rinse opened eye for several minutes under running water.
- **After swallowing:** Rinse out mouth and then drink plenty of water.
- **The following symptoms may occur:**  
After swallowing:  
Nausea  
Gastric pain  
Diarrhoea

(Contd. on page 2)

USA

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/03/2004

Reviewed on 02/03/2004

Trade name: *DPD No.1*

General feeling of sickness

(Contd. of page 1)

## 5 Fire fighting measures

- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **Special hazards caused by the material, its products of combustion or resulting gases:**  
Development of hazardous combustion gases or vapours possible in the event of fire.
- **Protective equipment:** Wear self-contained respiratory protective device.
- **Additional information:** Collect contaminated fire fighting water separately. It must not enter the sewage system.

## 6 Accidental release measures

- **Person-related safety precautions:** Not required.
- **Measures for environmental protection:** Do not allow product to reach sewage system or any water course.
- **Measures for cleaning/collecting:** Pick up mechanically.
- **Additional information:** No dangerous substances are released.

## 7 Handling and storage

- **Handling:**
- **Information for safe handling:** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a cool location.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Store in dry conditions.
- **Storage class:** Not required.

## 8 Exposure controls and personal protection

- **Additional information about design of technical systems:** No further data; see item 7.
- **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:** Preventive skin protection by use of skin-protecting agents is recommended.
- **Penetration time of glove material**  
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye protection:** Not required.
- **Body protection:** Protective work clothing

## 9 Physical and chemical properties

- **Odor Threshold:** Not applicable.
- **Form:** Tablets
- **Color:** White

(Contd. on page 3)

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/03/2004

Reviewed on 02/03/2004

Trade name: *DPD No.1*

(Contd. of page 2)

· <b>Odor:</b>	Odorless
· <b>Melting point/Melting range:</b>	Undetermined.
· <b>Boiling point/Boiling range:</b>	300°C (572°F)
· <b>Freezing Point:</b>	Not applicable.
· <b>Flash point:</b>	Not applicable.
· <b>Flammability (solid, gaseous):</b>	
<b>Upper Flammable Limit:</b>	Not applicable.
<b>Lower Flammable Limit:</b>	Not applicable.
· <b>Ignition temperature:</b>	Undetermined.
· <b>Sensitivity to Mechanical Impact:</b>	None
· <b>Sensitivity to Static Discharge:</b>	None
· <b>Auto igniting:</b>	Product is not selfigniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Vapor Pressure at 20°C (68°F):</b>	2.7 hPa (2 mm Hg)
· <b>Vapor Density:</b>	Not applicable.
· <b>Specific Gravity:</b>	Not applicable.
· <b>Density:</b>	Not determined.
· <b>Solubility in / Miscibility with Water:</b>	Soluble.
· <b>Coefficient of Water / Oil Distribution:</b>	Not applicable.
· <b>pH-value (10.5 g/l) at 20°C (68°F):</b>	6.3
· <b>Solvent content:</b>	
<b>Organic solvents:</b>	0.0 %
· <b>Solids content:</b>	100.0 %

## 10 Stability and reactivity

- **Thermal decomposition / conditions to be avoided:** strong heating
- **Dangerous reactions** No dangerous reactions known.
- **Dangerous products of decomposition:** No dangerous decomposition products known.

## 11 Toxicological information

### · Acute toxicity:

· <b>LD/LC50 values that are relevant for classification:</b>		
<b>10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight</b>		
Oral	LD50	2660 mg/kg (rat)
Dermal	LD50	> 2000 mg/kg (rat)
<b>7778-77-0 potassium dihydrogenorthophosphate</b>		
Oral	LD50	7100 mg/kg (rat)
<b>144-55-8 sodium hydrogencarbonate</b>		
Oral	LD50	4220 mg/kg (rat)
	LDLo	500 mg/kg (human)

- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.

(Contd. on page 4)

USA

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/03/2004

Reviewed on 02/03/2004

Trade name: **DPD No.1**

(Contd. of page 3)

· **Additional toxicological information:**

The product is not subject to classification according to internally approved calculation methods for preparations:  
When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

- **Carcinogenicity:** NTP? IARC Monographs? OSHA Regulated? Not found on these lists. see chapter 8 / 15
- **Teratogenicity:** Not found.
- **Mutagenicity:** Not found.
- **Reproductive Toxicity:** Not found.
- **Synergistic Products:** None

## 12 Ecological information

· **Ecotoxicological effects:**

· **Acquatic toxicity:**

**10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight**

EC50 | 133 mg/l/48h (Daphnia magna)

· **General notes:**

Water hazard class 2 (Self-assessment): hazardous for water  
Do not allow product to reach ground water, water course or sewage system.  
Danger to drinking water if even small quantities leak into the ground.  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.  
Must not reach bodies of water or drainage ditch undiluted or unneutralized.

## 13 Disposal considerations

· **Product:**

· **Recommendation:**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· **Uncleaned packagings:**

· **Recommendation:** Disposal must be made according to official regulations.

· **Recommended cleansing agent:** Water, if necessary with cleansing agents.

## 14 Transport information

· **DOT regulations:**

· **Hazard class:** None

· **Land transport ADR/RID (cross-border):**

· **ADR/RID class:** None

· **Maritime transport IMDG:**

· **IMDG Class:** None

· **Marine pollutant:** No

· **Air transport ICAO-TI and IATA-DGR:**

· **ICAO/IATA Class:** None

· **Propper shipping name:** None

· **TDG Class:** None

· **Transport/Additional information:** Not dangerous according to the above specifications.

USA

(Contd. on page 5)

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/03/2004

Reviewed on 02/03/2004

Trade name: **DPD No.1**

(Contd. of page 4)

## 15 Regulations

### · Sara

#### · Section 355 (Extremely hazardous substances):

None of the ingredients is listed.

#### · Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

#### · TSCA (Toxic Substances Control Act):

10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight

7778-77-0 potassium dihydrogenorthophosphate

144-55-8 sodium hydrogencarbonate

### · Proposition 65

#### · Chemicals known to cause cancer:

None of the ingredients is listed.

#### · Chemicals known to cause reproductive toxicity:

None of the ingredients is listed.

### · Canadian Ingredient Disclosure List

#### · Limit 0,1%

None of the ingredients is listed.

#### · Limit 1%

10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight

### · Other Regulations

#### · EPA (Environmental Protection Agency)

None of the ingredients is listed.

#### · IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

#### · NTP (National Toxicology Program)

None of the ingredients is listed.

#### · NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

#### · OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

### · Product related hazard informations:

Observe the general safety regulations when handling chemicals.

The product is not subject to identification regulations according to directives on hazardous materials.

#### · Water hazard class: Water hazard class 2 (Self-assessment): hazardous for water.

#### · CPR Classification: None

#### · MSDS Discloses Information Elements Required by the CPR .

## 16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

#### · Contact: Thermo Electron Corporation, EID, Water Analysis, Quality Assurance

#### · \* Data compared to the previous version altered.