

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/02/2004


Reviewed on 02/02/2004

## 1 Identification of substance

- **Product details**
- **Trade name:** Silica No.2
- **Product Use:** Reagent
- **Part of Orion:** AC2060
- **MSDS p/n:** 233266-001 Revision D
- **Manufacturer/Supplier:**  
Thermo Electron Corporation  
Environmental Instruments  
Water Analysis Group  
166 Cummings Center  
Beverly, MA 01915-6110  
USA  
Phone: 978-232-6000  
Made in Germany

## 2 Composition / Data on components

- **Chemical characterization**
- **Description:** Mixture: consisting of the following hazardous components.

CAS No.	Description	Hazard	Concentration
CAS: 87-69-4 EINECS: 201-766-0	(+)-tartaric acid	 Xi; R 36	40-45%

- **Additional information:**  
For the wording of the listed risk phrases refer to section 16.  
The remaining 55% contains no hazardous ingredients according to OSHA definition.

## 3 Hazards identification

- **Hazard description:**



Xi Irritant

- **Information pertaining to particular dangers for man and environment:**  
The product has to be labelled due to the calculation procedure of the "General Classification guideline for preparations of the EU" in the latest valid version.  
Irritating to eyes.
- **Classification system:**  
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.
- **NFPA ratings (scale 0 - 4)**  
Health = 1  
Fire = 1  
Reactivity = 0

## 4 First aid measures

- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Immediately rinse with water.
- **After eye contact:** Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**  
Drink copious amounts of water, induce vomiting.  
If symptoms persist consult doctor.

(Contd. on page 2)

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/02/2004

Reviewed on 02/02/2004

**Trade name: Silica No.2**

(Contd. of page 1)

**The following symptoms may occur:**

After swallowing:

- Nausea
- Diarrhoea
- Vomiting

**5 Fire fighting measures**

- **Suitable extinguishing agents:** Use fire fighting measures that suit the environment.
- **Special hazards caused by the material, its products of combustion or resulting gases:**  
Development of hazardous combustion gases or vapours possible in the event of fire.
- **Protective equipment:** Wear self-contained respiratory protective device.
- **Additional information** Collect contaminated fire fighting water separately. It must not enter the sewage system.

**6 Accidental release measures**

- **Person-related safety precautions:** Not required.
- **Measures for environmental protection:** Do not allow product to reach sewage system or any water course.
- **Measures for cleaning/collecting:** Pick up mechanically.

**7 Handling and storage**

- **Handling:**
- **Information for safe handling:** No special measures required.
- **Information about protection against explosions and fires:** No special measures required.
- **Storage:**
- **Requirements to be met by storerooms and receptacles:** Store in a cool location.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:**  
Keep receptacle tightly sealed.  
Store in dry conditions.
- **Storage class:** Not required.

**8 Exposure controls and personal protection**

- **Additional information about design of technical systems:** No further data; see item 7.
- **Components with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal protective equipment:**
- **General protective and hygienic measures:**  
Keep away from foodstuffs, beverages and feed.  
Immediately remove all soiled and contaminated clothing.  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes.  
Avoid contact with the eyes and skin.
- **Breathing equipment:** Use respiratory protective device against the effects of fumes/dust/aerosol.
- **Protection of hands:**  
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.  
After use of gloves apply skin-cleaning agents and skin cosmetics.

(Contd. on page 3)

USA

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/02/2004

Reviewed on 02/02/2004

Trade name: Silica No.2

(Contd. of page 2)

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Natural rubber, NR

PVC or PE gloves

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **Eye protection:** Tightly sealed goggles

· **Body protection:** Light weight protective clothing

## 9 Physical and chemical properties

· <b>Odor Threshold:</b>	Not applicable.
· <b>Form:</b>	Tablets
· <b>Color:</b>	Beige
· <b>Odor:</b>	Odorless
· <b>Melting point/Melting range:</b>	Undetermined.
· <b>Boiling point/Boiling range:</b>	300°C (572°F)
· <b>Freezing Point:</b>	Not applicable.
· <b>Flash point:</b>	> 100°C (> 212°F)
· <b>Flammability (solid, gaseous):</b>	
<b>Upper Flammable Limit:</b>	Not applicable.
<b>Lower Flammable Limit:</b>	Not applicable.
· <b>Ignition temperature:</b>	270.0°C (518°F)
· <b>Sensitivity to Mechanical Impact:</b>	None
· <b>Sensitivity to Static Discharge:</b>	None
· <b>Auto igniting:</b>	Product is not selfigniting.
· <b>Danger of explosion:</b>	Product does not present an explosion hazard.
· <b>Vapor Pressure at 20°C (68°F):</b>	2.7 hPa (2 mm Hg)
· <b>Vapor Density:</b>	Not applicable.
· <b>Specific Gravity:</b>	Not applicable.
· <b>Density at 20°C (68°F):</b>	1.656 g/cm <sup>3</sup>
· <b>Solubility in / Miscibility with</b>	
<b>Water:</b>	Soluble.
· <b>Coefficient of Water / Oil Distribution:</b>	Not applicable.
· <b>pH-value (11.6 g/l) at 20°C (68°F):</b>	3.2
· <b>Solvent content:</b>	
<b>Organic solvents:</b>	0.0 %
· <b>Solids content:</b>	100.0 %

## 10 Stability and reactivity

· **Thermal decomposition / conditions to be avoided:** strong heating

· **Materials to be avoided:** metals

· **Dangerous reactions** No dangerous reactions known.

· **Dangerous products of decomposition:** No dangerous decomposition products known.

USA

(Contd. on page 4)

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/02/2004

Reviewed on 02/02/2004

Trade name: Silica No.2

(Contd. of page 3)

## 11 Toxicological information

· **Acute toxicity:**

· **LD/LC50 values that are relevant for classification:**

**10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight**

Oral	LD50	2660 mg/kg (rat)
------	------	------------------

Dermal	LD50	> 2000 mg/kg (rat)
--------	------	--------------------

**134-03-2 sodium ascorbate**

Oral	LD50	11900 mg/kg (rat)
------	------	-------------------

· **Primary irritant effect:**

· **on the skin:** No irritant effect.

· **on the eye:** Irritating effect.

· **Sensitization:** No sensitizing effects known.

· **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:  
Irritant

· **Carcinogenicity:** NTP? IARC Monographs? OSHA Regulated? Not found on these lists. see chapter 8 / 15

· **Teratogenicity:** Not found.

· **Mutagenicity:** Not found.

· **Reproductive Toxicity:** Not found.

· **Synergistic Products:** None

## 12 Ecological information

· **Ecotoxicological effects:**

· **Aquatic toxicity:**

**10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight**

EC50	133 mg/l/48h (Daphnia magna)
------	------------------------------

· **Additional ecological information:**

· **BSB5-value:**

**87-69-4 (+)-tartaric acid**

BSB5	0.35 g/g (.)
------	--------------

· **General notes:**

Water hazard class 1 (Self-assessment): slightly hazardous for water

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

Must not reach bodies of water or drainage ditch undiluted or unneutralized.

## 13 Disposal considerations

· **Product:**

· **Recommendation:**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· **Uncleaned packagings:**

· **Recommendation:** Disposal must be made according to official regulations.

· **Recommended cleansing agent:** Water, if necessary with cleansing agents.

## 14 Transport information

· **DOT regulations:**

· **Hazard class:** None

(Contd. on page 5)

USA

# Material Safety Data Sheet

acc. to ISO/DIS 11014

Printing date 02/02/2004

Reviewed on 02/02/2004

Trade name: Silica No.2

(Contd. of page 4)

## · Land transport ADR/RID (cross-border):

· ADR/RID class: None

## · Maritime transport IMDG:

· IMDG Class: None

· Marine pollutant: No

## · Air transport ICAO-TI and IATA-DGR:

· ICAO/IATA Class: None

· Proper shipping name: None

· TDG Class: None

· Transport/Additional information: Not dangerous according to the above specifications.

## 15 Regulations

## · Sara

## · Section 355 (Extremely hazardous substances):

None of the ingredients is listed.

## · Section 313 (Specific toxic chemical listings):

None of the ingredients is listed.

## · TSCA (Toxic Substances Control Act):

87-69-4 (+)-tartaric acid

10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight

134-03-2 sodium ascorbate

## · Proposition 65

## · Chemicals known to cause cancer:

None of the ingredients is listed.

## · Chemicals known to cause reproductive toxicity:

None of the ingredients is listed.

## · Canadian Ingredient Disclosure List

## · Limit 0,1%

None of the ingredients is listed.

## · Limit 1%

10043-35-3 boric acid, crude natural, containing not more than 85 per cent of H3BO3 calculated on the dry weight

## · Other Regulations

## · EPA (Environmental Protection Agency)

None of the ingredients is listed.

## · IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

## · NTP (National Toxicology Program)

None of the ingredients is listed.

## · NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

## · OSHA-Ca (Occupational Safety &amp; Health Administration)

None of the ingredients is listed.

## · Product related hazard informations:

The product has been classified and marked in accordance with directives on hazardous materials.

(Contd. on page 6)

**Material Safety Data Sheet**  
**acc. to ISO/DIS 11014**

Printing date 02/02/2004

Reviewed on 02/02/2004

**Trade name: Silica No.2**

(Contd. of page 5)

**· Hazard symbols:***Xi Irritant***· Risk phrases:***Irritating to eyes.***· Safety phrases:***Keep out of the reach of children.**Avoid contact with eyes.**In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.**If swallowed, seek medical advice immediately and show this container or label.**Dispose of this material and its container at hazardous or special waste collection point**This material and its container must be disposed of as hazardous waste.***· Water hazard class:** *Water hazard class 1 (Self-assessment): slightly hazardous for water.***· CPR Classification:** *Class D, Div. 2***· MSDS Discloses Information Elements Required by the CPR .****16 Other information**

*This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.*

**· Contact:** *Thermo Electron Corporation, EID, Water Analysis, Quality Assurance***· \* Data compared to the previous version altered.**

—USA—