

# Certificate of Analysis

## Product Information:

**Orion Number** 910110  
**Description** # 10.01 BUFFER, 475 ML  
**Lot Number** EZ1  
**Expiration Date** 02/28/02  
**Method of Analysis** Potentiometric

## Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------|
| 1               | 10.018         | pH                     |

## Traceability:

| <u>N.I.S.T. Chemical</u>       | <u>SRM #</u> | <u>Sublot</u> |
|--------------------------------|--------------|---------------|
| Potassium Dihydrogen Phosphate | 186 I        | f             |
| Disodium Hydrogen Phosphate    | 186 II       | f             |
| Sodium Carbonate               | 192          | b             |
| Sodium Bicarbonate             | 191          | b             |

**Certificate Date:** 02/04/00

**Approved By:**

*Edward Henning*

**Quality Control Supervisor**

Environmental Instruments  
Water Analysis

166 Cummings Center

Beverly, MA  
01915

(978) 232-6015 Dom Fax  
(978) 232-6031 Int'l Fax

www.thermo.com  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

245967-001 Rev. B