

Certificate of Analysis

Product Information:

Orion Number 951211
Description ! AMMONIA PH ADJST ISA 475ML
Lot Number DQ1
Expiration Date 2 Years after Receipt of Product
Method of Analysis Potentiometric

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------------|
| 1 | 5.01 | Moles/Liter Sodium Hydroxide |
| 2 | 0.052 | Moles/Liter Disodium EDTA |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|------------------------------|--------------|---------------|
| Potassium Hydrogen Phthalate | 185 | g |

Certificate Date: 10/20/99

Approved By:

Edward Henning

Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245967-001 Rev. B

Thermo
ELECTRON CORPORATION