

# Certificate of Analysis

## Product Information:

**Orion Number** 910725  
**Description** PH 7 BUFFER IN 25 PACK.  
**Lot Number** CT2A  
**Expiration Date** 07/31/00  
**Method of Analysis** Potentiometric

## Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|------------------------|
| 1               | 7.002          | pH                     |

## Traceability:

| <u>N.I.S.T. Chemical</u>       | <u>SRM #</u> | <u>Sublot</u> |
|--------------------------------|--------------|---------------|
| Disodium Hydrogen Phosphate    | 186 II       | e             |
| Potassium Dihydrogen Phosphate | 186 I        | f             |
| Sodium Bicarbonate             | 191          | a             |
| Sodium Carbonate               | 192          | a             |

**Certificate Date:** 01/11/02

**Approved By:**



**Quality Control Supervisor**

Environmental Instruments  
Water Analysis

166 Cummings Center

Beverly, MA  
01915

(978) 232-6015 Dom Fax  
(978) 232-6031 Int'l Fax

www.thermo.com  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

245967-001 Rev. B