

Certificate of Analysis

Product Information:

Orion Number	967901
Description	ORP STANDARD, PINT BOTTLE
Lot Number	HZ1
Expiration Date	01/05
Method of Analysis	Potentiometric
Acceptance Range	7.050 +/- 0.015 pH Units/ORP Potential 220 +/- 1.0 mV/4.0 +/- 0.1 M Potassium Iodine/0.20 +/- 0.01 M Boric Acid/0.00069 +/- 0.0002 M Iodine

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	0.0069	Moles/Liter Iodine

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Dichromate	136	e

Certificate Date: 02/05/03

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

Thermo
ELECTRON CORPORATION

244721-001 Rev. B