

Certificate of Analysis

Product Information:

Orion Number	841113
Description	NA RECONDITIONING SOLN, 475 ML
Lot Number	HS1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	Fluoride 0.100 M +/- 0.005 M

Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	0.098	Moles/Liter Ammonium Bifluoride

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Fluoride	2203	NA

Certificate Date: 08/27/03

Approved By:

Edward Henning

Quality Control Supervisor