

Certificate of Analysis

Product Information:

Orion Number	967901
Description	ORP STANDARD, PINT BOTTLE
Lot Number	HS1
Expiration Date	8/05
Method of Analysis	Potentiometric
Acceptance Range	0.0069 +/- 0.0002 M Iodine

Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	0.0071	Moles/Liter Iodine

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Dichromate	136e	NA

Certificate Date: 09/02/03

Approved By:

Edward Henning

Quality Control Supervisor