

Certificate of Analysis

Product Information:

Orion Number	951207
Description	AMMONIA STD, 100PPM AS N 475ML
Lot Number	HS1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	100 ppm +/- 1 ppm Ammonia as Nitrogen

Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	99.5	ppm as Ammonium Chloride

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 08/29/03

Approved By:



Quality Control Supervisor