

Certificate of Analysis

Product Information:

Orion Number	951011
Description	! ALKALINE REAGENT 475 ML
Lot Number	IZ1
Expiration Date	2 Years after Receipt of Product
Method of Analysis	Potentiometric
Acceptance Range	10.00 N +/- 0.50 N Sodium hydroxide (NaOH)

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	9.95	Moles/Liter Sodium Hydroxide

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 01/12/04

Approved By:

Edward Henning

Quality Control Supervisor