

Certificate of Analysis

Product Information:

Orion Number	900022
Description	PH CLEANERS B, 4 X 60 ML
Lot Number	IZ1
Expiration Date	1/05
Method of Analysis	Potentiometric
Acceptance Range	2.5 % +/- 0.2 Bleach

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	2.5	% Weight by Volume Chlorine

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Dichromate	136e	NA

Certificate Date: 01/06/04

Approved By:

Edward Henning

Quality Control Supervisor