

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	910107
<b>Description</b>	# PH BUFFER 7, 475ML
<b>Lot Number</b>	HO3
<b>Expiration Date</b>	12/05
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	7.00 +/- 0.01 pH units

## Actual Analysis:


<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	7.009	pH @ 25 C

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	f
Potassium Dihydrogen Phosphate	186 I	f
Sodium Carbonate	192	b
Sodium Bicarbonate	191	b

**Certificate Date:** 12/31/03

**Approved By:**



**Quality Control Supervisor**