

Certificate of Analysis

Product Information:

Orion Number	700011
Description	TOTAL ALKALINITY REAGENT REFILL 4 X 475ML
Lot Number	IZ1
Expiration Date	1/05
Method of Analysis	Potentiometric
Acceptance Range	pH 3.50 +/- 0.05

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	3.51	pH @ 25 C

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	f
Potassium Dihydrogen Phosphate	186 I	f
Potassium Hydrogen Phthalate	185	g

Certificate Date: 01/14/04

Approved By:



Quality Control Supervisor